


Endobronchial Valves


A new solution for treating
persistent air leaks



Zephyr[®]
Endobronchial Valve

A Novel Approach to Treat Air Leaks

Data suggests air leaks may occur owing to the failure to fully seal lung tissue in up to 15% of lung surgery cases.¹ Leakage of air into the pleural space can lead to breathing difficulties and collapsed lungs.² Other causes of air leaks can include underlying pulmonary diseases, trauma, and spontaneous pneumothorax.³ Spontaneous pneumothorax occurs more frequently in patients with severe COPD.⁴ Other risk factors include diabetes, smoking, poor nutritional status and other underlying lung diseases.⁵



Once an air leak develops, it is possible to persist and worsen.³ Patients with persistent pneumothorax have significant comorbidities as shown by high American Society of Anesthesiologists scores.⁶ Complications of persistent pneumothorax can include pneumonia, empyema and fever.⁷ Patients with persistent pneumothorax therefore generally require extensive hospitalisation and resources, have a poor prognosis for surgery and have an urgent medical need for non-invasive techniques.^{1,2,9,10}

Zephyr to Block Air Flow

Once located, the air leak can be treated by implanting Endobronchial valve(s) (EBV(s)) using a delivery catheter.¹² The Zephyr EBVs block airflow to the area of leak in order to facilitate healing and allow the removal of chest tubes and patient discharge.^{3,12} Zephyr valves are easily removed using a bronchoscope.⁹




Clinical efficacy & safety

EBVs have been broadly utilized as a treatment option in this population resulting in a number of peer reviewed publications.^{1,8,9,11,13,14}

A series of 40 patients with persistent pneumothorax published in CHEST showed immediate improvement in 93 % of patients receiving Zephyr EBV.¹ Complete cessation of air leak was seen in 48 % of patients.

Another study looked at 16 patients with persistent pneumothorax (>7 days, various etiologies) who had not responded to standard treatment.⁹ Here again 77 % (10/13) of patients responded to valve therapy and mean air flow fell immediately from 871 ± 551 ml/min to 61 ± 72 ml/min when valves were implanted (see figure).




Mean air leak flow before and after valve implantation in a series of 16 patients non-responsive to standard treatment.


Mean implantation time was 30 ± 5 min to the end of the bronchoscopy.⁹ Valves were implanted in all lobes with equal efficacy.^{1,9}

EBVs are safe to use. In the two studies mentioned above, there were a total of 6 adverse events seen up to 42 months.^{1,9} None of the events was related to valve placement.

EBVs can remain in place, or can be removed without recurrence of pneumothorax. In the studies above, valves were removed from 15/56 patients, with no further complications.^{1,9}



Treatment algorithm¹



References

1. Travaline JM, McKenna RJ, De Giacomo T, Venuta F, Hazelrigg SR, Boomer M et al. Treatment of persistent pulmonary air leaks using endobronchial valves. *Chest* 2009;136(2):355-360.
2. Baumann MH. Pneumothorax. *Semin Respir Crit Care Med* 2001;22(6):647-656.
3. Miller S, Cohen J, Bellinger C, Chin R, Conforti J. Endobronchial Valves for the Treatment of Persistent Air Leaks. *Clin Pulm Med* 2012;19:143-146.
4. Rice TW, Okereke IC, Blackstone EH. Persistent air-leak following pulmonary resection. *Chest Surg Clin N Am* 2002;12(3):529-539.
5. Lois M, Noppen M. Bronchopleural fistulas: an overview of the problem with special focus on endoscopic management. *Chest* 2005;128(6):3955-3965.
6. Baumann MH, Strange C, Hefner JE, Light R, Kirby TJ, Klein J et al. Management of spontaneous pneumothorax: an American College of Chest Physicians Delphi consensus statement. *Chest* 2001;119(2):590-602.
7. Brunelli A, Xiume F, Al Refai M, Salati M, Marasco R, Sabbatini A. Air leaks after lobectomy increase the risk of empyema but not of cardiopulmonary complications: a case-matched analysis. *Chest* 2006;130(4):1150-1156.
8. Ambrosino N, Ribecchini A, Allidi F, Gabbrielli L. Use of endobronchial valves in persistent air leaks: a case report and review of the literature. *Expert Rev Respir Med* 2013;7(1):85-90.
9. Firlinger I, Stubenberger E, Muller MR, Burghuber OC, Valipour A. Endoscopic one-way valve implantation in patients with prolonged air leak and the use of digital air leak monitoring. *Ann Thorac Surg* 2013;95(4):1243-1249.
10. Cerfolio RJ, Bass CS, Pask AH, Katholi CR. Predictors and treatment of persistent air leaks. *Ann Thorac Surg* 2002;73(6):1727-31.
11. Ferguson JS, Sprenger K, Van Natta T. Closure of a bronchopleural fistula using bronchoscopic placement of an endobronchial valve designed for the treatment of emphysema. *Chest* 2006;129(2):479-481.
12. Anile M, Venuta F, De Giacomo T, Rendina EA, Diso D, Pugliese F et al. Treatment of persistent air leakage with endobronchial one-way valves. *J Thorac Cardiovasc Surg* 2006;132(3):711-712.
13. Brichon PY, Poquet C, Arvieux C, Pison C. Successful treatment of a life-threatening air leakage, complicating severe abdominal sepsis, with a one-way endobronchial valve. *Interact Cardiovasc Thorac Surg* 2012;15(4):779-780.
14. Weinreb N, Riker D, Beamis J, Lamb C. Ease of use of watanabe spigot for alveolopleural fistulas. *J Bronchology Interv Pulmonol* 2009;16(2):130-132.

