



Respiratory care accessories catalogue



Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Welcome to the ResMed respiratory care accessories catalogue

Thank you for choosing a ResMed ventilation solution. To help you maximise its potential, all our ventilators are available with a range of accessories, which can be used to further tailor therapy. Whether you are looking for added monitoring, transport options or humidification, this catalogue includes all the necessary product information to help you choose the right accessory for your applications. To keep updated with our latest available accessories, you can also visit ResMed.com*



S9 VPAP™

ResMed's VPAP series provides non-invasive bilevel devices for your non-dependent and spontaneously breathing patients with respiratory insufficiency. From initial setup to long-term care, the VPAP series delivers quality ventilation with a variety of accessories to help you meet a wide range of patient needs.



Lumis™

Lumis is a non-invasive ventilator series designed for non-dependent patients with respiratory insufficiency. Easy to set up and use, Lumis ventilators feature intelligent technologies, including iVAPS and AutoEPAP, which automatically adapt to individual breathing needs – making personalised ventilation easy.



Stellar™

The Stellar range is ResMed's premium non-invasive ventilator with invasive capabilities for non-dependent patients with mobility needs. This versatile device has a range of complimentary accessories, including a remote alarm, a mobility bag, internal and external battery options, and an SpO₂ and FiO₂ kit.



Astral™

Astral delivers excellence in leak and valve ventilation for invasive and non-invasive therapy. Featuring an impressive battery life, quick setup, and FiO₂ and SpO₂ monitoring options, it is suitable for a wide range of patients. It also offers mobility bags, external batteries and a dedicated hospital trolley for daily travel and convenience.












* Parts and accessories for Eliseé™ 250 and Eliseé 350 ventilators may still be available while stocks last. Please contact your local ResMed office or ResMed affiliate for more information.












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.














Contents
















Patient circuits, connectors and filters	2
Protective, transport and battery options	11
Power cords and clips	20
Monitoring tools and remote alarms	22
Data management	27
Humidifier accessories	28
Interfaces	30
Masks	33

Patient circuits, connectors, and filters		Product codes	 S9	 Lumis	 Stellar	 Astral
	SlimLine™ tubing (15 mm)	36810	•	•	•	
	ClimateLineAir™ heated tube	37296		•		
	ClimateLineAir Oxy tube	37357		•		
	ClimateLine™ tubing (15 mm)	36995	•			
	ClimateLine Max tubing (22 mm)	36997	•			
	ClimateLine Max Oxy tubing (22 mm)	36896 UK 36996	•			
	Standard tubing wrap	33963	•	•	•	

Patient circuits, connectors, and filters		Product codes	 S9	 Lumis	 Stellar	 Astral
	SlimLine tubing wrap	36811	•	•	•	
	Patient circuit, grey, autoclavable 22 mm x 2.0 m	14948	•	•	•	•
	Patient circuit, smooth bore: grey/clear with grey, non-autoclavable 22 mm x 2.0 m	21948	•	•	•	•
	Patient circuit, grey, autoclavable 22 mm x 3.0 m	14980	•	•	•	•
	Double limb circuit smooth bore with water trap 22 mm x 1.8 m	21301				•
	Double limb circuit smooth bore with water trap 15 mm x 1.6 m	21302				•
	Double limb circuit smooth bore 22 mm x 1.8 m Double limb circuit smooth bore 15 mm x 1.6 m	21303 21304				•












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Patient circuits, connectors, and filters		Product codes	 S9	 Lumis	 Stellar	 Astral
	Double limb circuit smooth bore 10 mm x 1.6 m	21305				•
	Double limb circuit corrugated with water trap 22 mm x 1.6 m Double limb circuit corrugated with water trap 15 mm x 1.6 m	21306 21307				•
	Double limb circuit corrugated 22 mm x 1.6 m Double limb circuit corrugated 15 mm x 1.6 m	21308 21309				•
	Double limb circuit coaxial 30 mm x 1.6 m	21323				•
	Single limb circuit smooth bore with exh valve 22 mm x 1.8 m* Single limb circuit smooth bore with exh valve 15 mm x 1.6 m*	21311 21312				•
	Single limb circuit smooth bore with water trap and exh valve 22 mm x 1.8 m* Single limb circuit smooth bore with water trap and exh valve 15 mm x 1.6 m*	21313 21314				•
	Single limb circuit corrugated with exh valve 22 mm x 1.6 m* Single limb circuit corrugated with exh valve 15 mm x 1.6 m*	21321 21322				•












Patient circuits, connectors, and filters	Product codes	 S9	 Lumis	 Stellar	 Astral
 <p>Single limb circuit corrugated with water trap and exh valve 22 mm x 1.6 m*</p>  <p>Single limb circuit corrugated with water trap and exh valve 15 mm x 1.6 m*</p>	<p>21324 21325</p>				•
 <p>Double limb corrugated heated 22 mm x 1.6 m</p>  <p>Double limb corrugated heated 15 mm x 1.6 m</p>	<p>21332 21334</p>				•
 <p>Single limb corrugated heated with exh valve 22 mm x 1.6 m*</p>  <p>Single limb corrugated heated with exh valve 15 mm x 1.6 m*</p>	<p>21335 21336</p>				•
 <p>Catheter mount smooth bore with port 22F</p>	<p>21341</p>				•
 <p>Catheter mount corrugated with port 22F</p>	<p>21342</p>				•
 <p>Catheter mount paediatric 15M</p>	<p>21343</p>				•
 <p>Double limb corrugated heated with F&P chamber 22 mm x 1.6 m</p>  <p>Double limb corrugated heated with F&P chamber 15 mm x 1.6 m</p>	<p>21331 21333</p>				•

*These single limb circuits include ResMed's QuickConnect™ with expiratory valve.

Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.












Patient circuits, connectors, and filters		Product codes	 S9	 Lumis	 Stellar	 Astral
	Double limb adapter	27940				•
	Double limb adapter – Clear Cap	27980				•
	Single limb valve adapter*	27941				•
	Single limb leak adapter*	27942				•
	Air10™ tubing elbow	37394		•		
	Air10 outlet	37310		•		
	ResMed leak valve Note: this replaces ResMed leak port 24976	24988			•	•












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Patient circuits, connectors, and filters		Product codes	 S9	 Lumis	 Stellar	 Astral
	Oxygen connector port (100 pack)	26960	•	•		
	Anti-bacterial filter (1 pack) Anti-bacterial filter (50 pack)	24966-1 24966	•	•	•	•
	Heat and moisture exchange filter (HMEF) for adults (1 pack) HMEF for adults (50 pack)	24967-1 24967			•	•
	HMEF for paediatric patients (1 pack) HMEF for paediatric patients (50 pack)	24969-1 24969			•	•
	HME maxi for adults (50 pack) HME midi for paediatric patients (50 pack)	27920 27921				•
	Stellar filter cover (ROW)	24885			•	
	VPAP filter standard (1 pack) VPAP filter standard (2 pack) VPAP filter standard (12 pack) VPAP filter standard (50 pack)	36850 36851 36852 36853	•	•		











*These single limb circuits include ResMed's QuickConnect™ with expiratory valve.











Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Patient circuits, connectors, and filters		Product codes	 S9	 Lumis	 Stellar	 Astral
	VPAP filter hypo-allergenic (1 pack) VPAP filter hypo-allergenic (2 pack) VPAP filter hypo-allergenic (12 pack) VPAP filter hypo-allergenic (50 pack)	36855 36856 36857 36858	•	•		
	VPAP filter cover	36859	•			
	Stellar filter hypo-allergenic (1 pack) Stellar Filter hypo-allergenic (12 pack)	24933 24934			•	
	Air10 filter cover – Charcoal	19729		•		
	Astral inlet filter (4 pack)	27939				•
	HMEF for paediatric patients – small angled (50 pack)	24968			•	•
	Male-male connector 22 mm	707322	•	•	•	•











Protective, transport, and battery options		Product codes	 S9	 Lumis	 Stellar	 Astral
	ResMed Power Station II carry bag	24952	•	•	•	
	VPAP carry bag	36862	•			
	Air10 travel bag	37304		•		
	Stellar carry bag	24989			•	
	Stellar mobility bag	24977			•	
	Astral mobility bag	27917				•
	SlimFit mobility bag	27949				•









Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Protective, transport, and battery options		Product codes	 S9	 Lumis	 Stellar	 Astral
	Hospital trolley (includes base plate)	27914				•
	Hospital trolley storage bin	27929				•
	Hospital trolley cylinder mount	27930				•
	Hospital trolley power strip - EU Hospital trolley power strip - ANZ Hospital trolley power strip - UK	27932 27933 27934				•
	ResMed standard trolley (does not include top plate)	24905	•	•	•	
	Circuit support arm, compatible with standard trolley 24905 Tube clip, single, for circuit support arm 24917 Tube clip, double for circuit support arm 24917	24917 24920 24919	•	•	•	












Protective, transport, and battery options		Product codes	 S9	 Lumis	 Stellar	 Astral
	Astral table stand	27915				•
	Astral hospital trolley base plate	27928				•
	Top plate for standard trolley	24913	•	•	•	
	Astral external battery	27918				•
	ResMed Power Station II Europe and APAC – includes carry bag and power supply adaptor (does not include DC cable or power supply). EU: 36823 UK: ANZ: 36824	24921	•	•	•	•
	ResMed Power Station II power supply unit (90W) including power cord: (For non-S9 series users.)					
	Power cord Japan	36825		•	•	
	Power cord UK	36822				
	Power cord Europe	36823				
	Power cord ANZ	36824				












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Protective, transport, and battery options		Product codes	 S9	 Lumis	 Stellar	 Astral
	ResMed Power Station II power supply unit adapter	24946	•	•	•	•
	ResMed Power Station II	24926			•	•
	ResMed Power Station II DC cable (0.5 m) for use in the Stellar mobility bag (For non-S9 series users.)	24981			•	
	ResMed Power Station II DC cable (For non-S9 series users)	24959			•	
	ResMed Power Station II coupler kit to connect two batteries	24964	•	•	•	
	MPV Circuit Support Arm	27955	•	•	•	•












Protective, transport, and battery options		Product codes	 S9	 Lumis	 Stellar	 Astral
	Homecare Stand – Astral (bundle)					
	Homecare Stand – S9 with Humidifier (bundle)	279001				
	Homecare Stand - S9 (bundle)	279002				
	Homecare Stand - Air 10 (bundle)	279003	•	•		•
	Homecare Stand – Lumis TX kit (bundle)	279004				
	Homecare Stand – Lumis TX Humidifier kit (bundle)	360121 360122				
Note: Devices not included in the bundle.						
	Homecare Stand without Interface/Top Plate	27963				
	Astral Top Plate	27956				
	S9 Series with Humidifier Top Plate	27957				
	S9 Series without Humidifier Top Plate	27958	•	•		•
	Air 10 – Lumis Series Top Plate (ROW)	27964				
	Homecare Stand Screws spare kit	27960				
	Homecare Stand Casters Spare Part	27962				
Note: Representative image shown.						
	Homecare Stand Astral Plate (attaches to bottom of Astral)	27961				•
	Homecare Stand IV Pole Kit	27948				











Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Power cords and clips		Product codes	 S9	 Lumis	 Stellar	 Astral
	Air10 90W PSU Air10 90W PSU UK Air10 90W PSU Europe	37347 37346 37345		•		
	Air10 DC/DC converter	37297		•		
	VPAP DC/DC converter (24V/90W)	36970	•			
	Stellar DC/DC converter (24V/90W)	24922			•	
	Power cord clip (2 pack)	26914			•	
	Stellar power cord Japan Stellar power cord APAC Stellar power cord Europe Stellar power cord UK Stellar power cord Korea	24998 24973 24974 24975 24997			•	
	Power supply unit (90W) with power cord Japan Power supply unit (90W) with power cord UK Power supply unit (90W) with power cord Europe Power supply unit (90W) with power cord ANZ Power supply unit (30W) without power cord Power supply unit (30W) without power cord Japan	36825 36822 36823 36824 36920 36921			•	•








Monitoring tools and remote alarms		Product codes	 S9	 Lumis	 Stellar	 Astral
	Remote Alarm II	27902			•	•
	Remote Alarm II cable (2 m) Remote Alarm II cable (5 m) Remote Alarm II cable (10 m) Remote Alarm II cable (20 m)	27904 27905 27906 27907			•	•
	Hospital alarm system cable (30 m)	27912			•	•
	Oxygen monitoring kit (does NOT include oxygen sensor CAP008617) with connecting FiO ₂ cable (CAB015770) and T piece adapter (RAC008157) for single patient use	24955			•	
	Astral oxygen sensor	27935				•
	Astral expiratory flow sensor	27936				•
	Low pressure oxygen connector	24971			•	•












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Monitoring tools and remote alarms	Product codes	 S9	 Lumis	 Stellar	 Astral	
	FiO ₂ cable (for oxygen sensor CAP008617)	CAB015770			•	
	T-piece adapter for FiO ₂ monitoring kit	RAC008157			•	
	FiO ₂ monitoring sensor (CAP008617)	24980			•	
	Nonin soft sensor (medium)	70568		•	•	•
	Oximetry module	37302		•		
	VPAP oximeter adapter	36940	•			
	Oximeter flex sensor (25 pack) (3 m, 1 sensor, 25 tapes)	16936	•	•	•	•










Monitoring tools and remote alarms		Product codes	 S9	 Lumis	 Stellar	 Astral
	VPAP oximeter converter cable	36941	•			
	VPAP oximeter extension cable. The extension cable connects to the sensor on one end, and the XPOD on the other end	1430500	•			
	Oximetry accessories pack (RoHS compliant) includes: ResMed Xpod oximeter (22305), clip for Xpod (22306) and medium soft sensor 3 m cable (707560) Note: this pack replaces product code 36942	36944	•	•		
	ResMed Xpod 3012 LP oximeter (replaces Xpod oximeter 22305) Nonin Xpod clip Oximetry set with reusable soft sensor (includes Xpod LP oximeter and clip plus a Nonin 8000SM reusable soft sensor) Oximetry set with disposable sensor (includes Xpod LP oximeter and clip plus 3 x Nonin 7000A single-use oximetry sensors)	22374 22371 22370 22380	•	•	•	•
	ResMed Connectivity Module (RCM)	27201 (US) 27202 (EUR) 27203 (APAC)			•	•
	RCM USB Cable	27953			•	•

Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Monitoring tools and remote alarms		Product codes	 S9	 Lumis	 Stellar	 Astral
	RCM Power Supply	27954			•	•
	RCM Wall Mount	27975			•	•
	RCM Stand	27976			•	•












Data management		Product codes	 S9	 Lumis	 Stellar	 Astral
	Embletta adapter for the VPAP series	36955	•			
	EasyView software MKG014454 and serial line cord CAB012527	KIT014618				
	ResScan data management software for data download. Visit ResMed.com/resscan/compatibility for data management compatibility	31302	•	•	•	•
	USB cable (connecting the device to a PC)	629051	•	•	•	•
	ResMed USB	24935			•	•
	SD card with envelopes Eur Group I (1 pack) SD card with envelopes Eur Group II (1 pack) SD card with envelopes Eur Group I (10 pack) SD card with envelopes Eur Group II (10 pack)	36843 36844 36841 36842	•			
	SD card reader for VPAP series	36931	•			

Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.












Data management		Product codes	 S9	 Lumis	 Stellar	 Astral
	USB adapter for VPAP series	36950	•			
	USB extension cable for VPAP series	7071015	•			
	NIV+ with user guides and stickers (10 pack)**	KIT015934				•
	USB module	37301		•		
	SD card with protective folder – (1 pk)	37334 37359		•		
	SD card with protective folder – (10 pk)	37332 37358				







**English, Polish, Russian, Czech, Slovenian, Hungarian








Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Humidifier accessories		Product codes	 S9	 Lumis	 Stellar	 Astral
	H4i™ integrated heated humidifier Germany H4i integrated heated humidifier Eur Group I and UK H4i integrated heated humidifier Eur Group II and APAC H4i integrated heated humidifier ANZ H4i integrated heated humidifier France	24927 24929 24930 26951 24936			•	
	Stellar H4i water chamber (country dependent, disposable)	26952			•	
	Stellar H4i cleanable water chamber with spill guard	26958			•	
	Stellar H4i top cover and seal	26953			•	
	Stellar H4i cleaning plug	26955			•	
	Stellar H4i tub plate seal	26954			•	
	Stellar H4i connector plug (10 pack)	24987			•	












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Humidifier accessories		Product codes	 S9	 Lumis	 Stellar	 Astral
	H5i™ integrated heated humidifier Japan H5i integrated heated humidifier Eur Group I H5i integrated heated humidifier ROW H5i integrated heated humidifier France H5i integrated heated humidifier Germany	36900 36904 36905 36906 36907	•			
	H5i cleanable water chamber for VPAP series	36800	•			
	H5i water chamber for VPAP series	36803	•			
	H5i flip lid for VPAP series	36891	•			
	H5i flip lid seal for VPAP series	36892	•			
	H5i outlet clip for the VPAP series	36894	•			
	HumidAir™ heated humidifier – Cleanable tub (comes with device as standard)	37300		•		

Humidifier accessories		Product codes	 S9	 Lumis	 Stellar	 Astral
	HumidAir Standard Tub	37299		•		
	Air10 Side Cover Charcoal	37303		•		












Interfaces		Product codes	 S9	 Lumis	 Stellar	 Astral
	EasySpeak™ mouthpiece angled for 15 M connection	21351				•
	EasySpeak mouthpiece angled for 22 M/15 F connection	21354				•
	EasySpeak mouthpiece angled with 17 cm flexi-tube 22 F	21353				•

Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.












Masks*	Mask types	 S9	 Lumis	 Stellar	 Astral
 AirFit™ F20	Full Face Mask	•	•	•	•
 AirFit™ F10	Full Face Mask	•	•	•	•
 Quattro™ Air (validated up to 40 cm)	Full Face Mask	•	•	•	•
 Quattro Air For Her (validated up to 40 cm)	Full Face Mask	•	•	•	•
 Quattro FX	Full Face Mask	•	•	•	•
 Quattro FX for Her	Full Face Mask	•	•	•	•
 Mirage Quattro™ (validated up to 40 cm)	Full Face Mask	•	•	•	•

* Please refer to ResMed.com/masks for more details on masks












Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Masks*	Mask types	 S9	 Lumis	 Stellar	 Astral
 Ultra Mirage™	Full Face Mask	•	•	•	•
 Mirage Liberty™	Full Face Mask	•	•	•	•
 Mirage Micro™	Nasal Mask	•	•	•	•
 Mirage Activa™	Nasal Mask	•	•	•	
 Mirage Activa LT	Nasal Mask	•	•	•	
 Ultra Mirage™ II	Nasal Mask	•	•	•	•
 Pixi™	Nasal Mask	•	•	•	•











Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.







Masks*		Mask types	 S9	 Lumis	 Stellar	 Astral
	Mirage SoftGel	Nasal Mask	•	•	•	•
	Mirage FX	Nasal Mask	•	•	•	•
	Mirage FX for Her	Nasal Mask	•	•	•	•
	AirFit™ N20	Nasal Mask	•	•	•	•
	AirFit N10	Nasal Mask	•	•	•	•
	Swift™ FX Nano	Nasal Mask	•	•	•	•
	Swift FX Nano for Her	Nasal Mask	•	•	•	•

Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Masks*	Mask types	 S9	 Lumis	 Stellar	 Astral
 AirFit P10	Nasal Pillows Mask	•	•	•	
 Swift FX Bella	Nasal Pillows Mask	•	•	•	•
 Swift FX Bella Gray	Nasal Pillows Mask	•	•	•	•
 Mirage Swift™ II	Nasal Pillows Mask	•	•	•	•
 Swift LT	Nasal Pillows	•	•	•	•
 Swift FX	Nasal Pillows	•	•	•	•
 Swift FX for Her	Nasal Pillows	•	•	•	•


Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.

Masks*	Mask types	 S9	 Lumis	 Stellar	 Astral
 Ultra Mirage NV full face	Non-Vented Full Face Mask				•
 Quattro FX NV	Non-Vented Full Face Mask				•
 Ultra Mirage NV nasal mask	Non-Vented Nasal Mask				•
 AcuCare™ F1-0	Hospital Full Face Mask Non-Vented				•
 AcuCare F1-1	Hospital Full Face Mask Vented with Anti-Asphixia Valve (AAV)				•
 AcuCare F1-4	Hospital Full Face Mask Vented	•	•	•	•

Masks*	Mask types	 S9	 Lumis	 Stellar	 Astral
 AcuCare HFNC	High Flow Nasal Cannula (HFNC)				
 ResMed Hospital Nasal	Hospital Nasal Mask			•	•

Note: Some products and accessories may not have a complete list of compatible ResMed ventilators indicated within this catalogue due to ongoing compatibility testing. Please consult the individual product user guides and compatibility guides for conditions of use, warnings and precautions associated with each listed item. Photos included in this catalogue may not be an exact representation of the final product and not all products or accessories are available in all regions.



 **ResMed Ltd.**
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia
MANUFACTURER

CE 0123

Distributed by: ResMed Ltd, Bella Vista, NSW, Australia. See ResMed.com for ResMed locations worldwide. Activa, Astral, ClimateLine, ClimateLineAir, EasySpeak, Elisée, H5i, H4i, HumidAir, Mirage, Mirage Activa, Mirage Liberty, Mirage Micro, Mirage Quattro, Mirage Swift, Pixi, Quattro, ResScan, SiimLine, Stellar, Swift, TiControl, Ultra Mirage and VPAP are trademarks and/or registered trademarks of the ResMed family of companies. Nonin and Xpod are trademarks of Nonin Medical. For patent information, see ResMed.com/ip. Product availability may vary across regions. Specifications may change without notice. © 2017 ResMed Ltd. 1013296/13 2017-07

ResMed.com