

Seeing surgery through  
new eyes...

## IDEAL EYES

- Enhanced center-to-edge resolution
- Improved depth of field
- Reduced image distortion



# Seeing surgery through new eyes...

## IDEAL EYES

### The Way Surgery Was Meant to be Seen

Stryker is one of the visualization leader in high definition imaging with improved detail recognition, resolution, and excellent durability.

Designed to enhance picture quality, IDEAL EYES scopes are designed to:

- Enhance picture quality
- Enhance resolution
- Enhance colour reproduction
- Complement Stryker's latest camera technology: 1288 HD camera system

### Take A Better Look

IDEAL EYES scopes have an excellent optical performance because:

- 1) Adjusted optical calculations
- 2) High quality glass
- 3) US Patented lens design

### Set Your Sights High

Stryker IDEAL EYES scopes combine high quality sapphires and laser welded leakage paths to ensure durability that will withstand repeated autoclave cycles over time.



Stryker's latest HD endoscope generation includes:

- Laparoscopes 2.9mm, 5.5mm and 10mm
- Cystoscopes 4mm
- Hysteroscopes 2.9mm and 4mm
- Arthroscopes 1.9mm, 2.3mm, 2.7mm and 4mm

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Ideal Eyes, Stryker. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTXENDIEBREN  
MTX5713/GS 02/10

Copyright © 2010 Stryker



Stryker SA  
Cité Centre  
Grand-Rue 90  
1820 Montreux  
Switzerland

t: +41 21 966 12 01  
f: +41 21 966 12 00

[www.stryker.eu](http://www.stryker.eu)