

Excellent color reproduction
Greater clarity
Sharp images
Improved depth perception
Ergonomic camera head

PROCAM HD HIGH DEFINITION CAMERA

Ecleris continues to advance, introducing the newest addition to our endoscopic product line, ProCam HD



PROCAM HD HIGH DEFINITION CAMERA

PROCAM HD HIGH DEFINITION CAMERA



The Ecleris ProCam HD is a single chip CCD with 1280 x 720 lines of resolution and produces clear sharp images that were once only available on higher priced systems found in the surgical suite. Advanced features on the ProCam HD include Anti-moire (for flexible fiber optic endoscopes) and stroboscopy (a feature which balances the light used when performing video stroboscopy). With 3 HDMI outputs you can easily run multi video monitor configurations from a single platform.



- On screen display
 - White Balance
 - Windows function (2 levels)
 - Anti – moiré (2 levels)
 - Stroboscopy setup
 - Boost (Gain control)

TECHNICAL SPECIFICATIONS

Resolution	1280 X 720 p
Video System	progressive scan 60 / 50hz
CCD Sensor	1/3" HD
Sensitivity	5 LUX
Settings	on screen display
Outputs	3 X HDMI 1 X COMPOSITE (NTSC/PAL)
Power supply	100 – 240 VAC
Dimensions (w x d x h)	300 X 290 X 95 cm.
Head	C – mount thread
Class	I
Type	B

Rev. 10.08.12



www.ecleris.com

ECLERIS
INNOVATION VALUE