

REF. #950-001

DragonFly™ Surgical Drill System

The system includes the following:

ITEM	REF. #
Handpiece	950-002
Footswitch	950-003
Rechargeable Batteries	950-004
Battery Charger	950-005
Angled Adaptor, 15°	950-006-15
Angled Adaptor, 25°	950-006-25
Angled Adaptor, 35°	950-006-35
Sterilization Tray	950-007

Disposable burs available on demand:

ITEM	REF. #
0.4mm, 95mm L, Diamond	951-95-004
0.5mm, 95mm L, Diamond	951-95-005
0.6mm, 95mm L, Diamond	951-95-006
0.7mm, 95mm L, Diamond	951-95-007
0.8mm, 95mm L, Diamond	951-95-008
1.0mm, 95mm L, Diamond	951-95-010
1.0mm, 95mm L, Diamond w/ BurShield™	951-95-010-SL
1.4mm, 95mm L, Diamond	951-95-014
1.4mm, 95mm L, Diamond w/ BurShield™	951-95-014-SL
1.8mm, 95mm L, Diamond	951-95-018
1.8mm, 95mm L, Diamond w/ BurShield™	951-95-018-SL
2.3mm, 95mm L, Diamond	951-95-023
2.3mm, 95mm L, Diamond w/ BurShield™	951-95-023-SL
0.4mm, 80mm L, Carbide w/ BurShield™	952-80-004-SL
0.5mm, 80mm L, Carbide w/ BurShield™	952-80-005-SL
0.6mm, 80mm L, Carbide w/ BurShield™	952-80-006-SL
0.7mm, 80mm L, Carbide w/ BurShield™	952-80-007-SL
0.8mm, 80mm L, Carbide w/ BurShield™	952-80-008-SL
1.0mm, 80mm L, Carbide w/ BurShield™	952-80-010-SL
1.4mm, 80mm L, Carbide w/ BurShield™	952-80-014-SL
1.8mm, 80mm L, Carbide w/ BurShield™	952-80-018-SL
2.3mm, 80mm L, Carbide w/ BurShield™	952-80-023-SL
1.8mm, 80mm L, Egg w/ BurShield™	953-80-018-SL
2.3mm, 80mm L, Egg w/ BurShield™	953-80-023-SL
1.8mm, 80mm L, Cylinder w/ BurShield™	954-80-018-SL
0.35mm, 95mm L, Trocar Perforator	955-95-0035



Grace Medical We Listen.


8500 Wolf Lake Drive, Suite 110
Memphis, TN 38133 USA
Phone 1.901.386.0990
Fax 1.901.386.0950
Toll Free 1.866.472.2363
info@gracemedical.com
www.gracemedical.com
LIT0172 CID6035 Rev. 2018-09



DragonFlyTM

Surgical Drill System

powered by **InvisiDriveTM**



with patented
BurShieldTM technology



Grace Medical We Listen.



DragonFly™

Surgical Drill System




BurShield™

Adjustable shield for protection of nerve, soft tissue and intraoperative dressing




Grace Medical We Listen.



InvisiDrive™

Magnetic coupling for quick and easy exchange of burs



Angled Adaptors

Adaptors to vary bur angulation, maximizing visualization (15°, 25°, 35°)



Powered by InvisiDrive™, the DragonFly™ Surgical Drill System represents a model of user-friendly ergonomics. The handpiece is small in size, light in weight, and hermetically sealed to ensure optimal and consistent performance.

