


SEPLOU[®]
MEDICAL



Urology

Practical | Professional | Minimally Invasive

SEPLOU, Inc.

Since the company's founding in 2014, SEPLOU has been an innovator in the urological surgical medical technology. This is supported by a core group of experienced business leaders and medical professionals and physicians.

SEPLOU supports the work of medical professionals by providing advanced, minimally invasive therapeutic and diagnostic technologies to improve the quality of patient care around the globe.

Table of Contents

Ureteral Stents Set – Standard Double –J	01-02
Ureteral Stents Set –Standard Single-J	03
Ureteral Access Sheath	04
PTFE Guidewire	05
Hydrophilic Guidewire	05
Zebra Guidewire	06
Percutaneous Nephrostomy Set	07-08
Ureteral Catheter	09
Puncture Needles	10
Endoscope Protection Cover	11
Surgical Adhesive Membrane	11
Forceps	12
Electrodes	13-14
Silicone Foley Catheter 2-Way/3-Way	15
Latex Foley Catheter 2-Way/3-Way	16
Urine Bag	17
Leg Bag	18
Intermittent Catheter	19-20
Hydrophilic Intermittent Catheter	20-21

Ureteral Stents set

- Standard Double-J
- Multiloop

Made of Polyurethane



Ureteral Stents set

Set consists of stent, pusher

Standard Double-J

Closed end	Open end	Size(Fr)	Length(cm)
/	SPUS2540	4Fr	16
SPUS2140	SPUS2340	4Fr	26
SPUS2147	SPUS2347	4.7Fr	26
SPUS2160	SPUS2360	6Fr	26
SPUS2170	SPUS2370	7Fr	26
SPUS2180	SPUS2380	8Fr	26
SPUS2240	SPUS2440	4Fr	28
SPUS2247	SPUS2447	4.7Fr	28
SPUS2260	SPUS2460	6Fr	28
SPUS2270	SPUS2470	7Fr	28
SPUS2280	SPUS2480	8Fr	28

Multiloop

Closed end	Open end	Size(Fr)	Length(cm)
SPUS2740	SPUS2840	4Fr	Multiloop
SPUS2747	SPUS2847	4.7Fr	Multiloop
SPUS2760	SPUS2860	6.0Fr	Multiloop
SPUS2770	SPUS2870	7.0Fr	Multiloop
SPUS2780	SPUS2880	8.0Fr	Multiloop

Standard Double-J with PTFE Guidewire

Open end	Size(Fr)	Length(cm)	Guidewire
SPUS3340	4Fr	26	ø0.028"
SPUS3347	4.7Fr	26	ø0.035"
SPUS3360	6.0Fr	26	ø0.035"
SPUS3370	7.0Fr	26	ø0.035"
SPUS3380	8.0Fr	26	ø0.035"
SPUS3440	4Fr	28	ø0.028"
SPUS3447	4.7Fr	28	ø0.035"
SPUS3460	6.0Fr	28	ø0.035"
SPUS3470	7.0Fr	28	ø0.035"
SPUS3480	8.0Fr	28	ø0.035"

Ureteral Stents Set

-Single J

Made of polyurethane



Standard Single-J

Closed end	Open end	Size(Fr)	Length(cm)
SPUS1970	SPUS11070	7	50
SPUS1980	SPUS11080	8	50

Standard Single-J with PTFE Guidewire

Closed end	Open end	Size(Fr)	Length(cm)	Guidewire
SPUS4960	SPUS41060	6	50	ø0.035"
SPUS4970	SPUS41070	7	50	ø0.035"
SPUS4980	SPUS41080	8	50	ø0.035"

Ureteral Access Sheath

Provides ureteral dilation and a continuous working channel for the introduction of Endoscopes and Instruments during ureteral access procedures.



Product Codes	Length(cm)	Size (Fr)
SPAS1035	35	10Fr
SPAS1135	35	11Fr
SPAS1235	35	12Fr
SPAS1335	35	13Fr
SPAS1435	35	14Fr
SPAS1535	35	15Fr
SPAS1045	45	10Fr
SPAS1145	45	11Fr
SPAS1245	45	12Fr
SPAS1345	45	13Fr
SPAS1445	45	14Fr
SPAS1545	45	15Fr

PTFE Guidewire


PTFE coated

one tip stiff - one tip flexible

Colour: PTFE dark green

Made of stainless steel

Product Codes	Length(cm)	Size	Tip
SPGW1118	150	ø0.018"	J
SPGW1125	150	ø0.025"	J
SPGW1128	150	ø0.028"	J
SPGW1132	150	ø0.032"	J
SPGW1135	150	ø0.035"	J
SPGW1138	150	ø0.038"	J
SPGW1218	150	ø0.018"	Straight
SPGW1225	150	ø0.025"	Straight
SPGW1228	150	ø0.028"	Straight
SPGW1232	150	ø0.032"	Straight
SPGW1235	150	ø0.035"	Straight
SPGW1238	150	ø0.038"	Straight



Hydrophilic Guidewire

Hydrophilic coated

one tip stiff - one tip flexible

Colour: black

Made of Nitinol

Product Codes	Length(cm)	Size	Tip
SPGW2132	150	ø0.032"	Angled
SPGW2135	150	ø0.035"	Angled
SPGW2138	150	ø0.038"	Angled
SPGW2232	150	ø0.032"	Straight
SPGW2235	150	ø0.035"	Straight
SPGW2238	150	ø0.038"	Straight




Zebra Guidewire

65cm Hydrophilic coated

one tip stiff - one tip flexible

Made of Nitinol



Product Codes	Length(cm)	Size	Colour	Tip
SPGW3132	150	ø0.032"	Black & White	Straight
SPGW3135	150	ø0.035"	Black & White	Straight
SPGW3332	150	ø0.032"	Black & Yellow	Hydrophilic coated
SPGW3335	150	ø0.035"	Black & Yellow	Hydrophilic coated

Percutaneous Nephrostomy Set



Percutaneous Nephrostomy Set

Features:

- Hydrophilic coated sheath reduces patient's trauma
- Ease of Access promoted by gradual taper tip
- Excellent Radiopacity for improved visualization
- Available in 8 Fr to 20 Fr sizes in standard lengths

Includes:

- Dilator with Sheath (size and configuration may vary)
- Dilator (size and configuration may vary)
- Trocar (18G/20G)
- * Guidewire (optional)
- * 100% Silicone Integrat Catheter (optional)

Percutaneous Nephrostomy Set

Product Codes	Sheath Size (Fr)	Dilators Size (Fr)
SPPN0208	8Fr	8Fr
SPPN0210	10Fr	10Fr
SPPN0212	12Fr	8Fr-12Fr
SPPN0214	14Fr	8Fr-14Fr
SPPN0216	16Fr	8Fr-16Fr
SPPN0218	18Fr	8Fr-18Fr
SPPN0220	20Fr	8Fr-20Fr
SPPN0222	22Fr	8Fr-22Fr

Set with PTFE Guidewire

Product Codes	Sheath Size (Fr)	Dilators Size (Fr)	Guidewire
SPPN0308	8Fr	8Fr	ø0.035"
SPPN0310	10Fr	10Fr	ø0.035"
SPPN0312	12Fr	8Fr-12Fr	ø0.035"
SPPN0314	14Fr	8Fr-14Fr	ø0.035"
SPPN0316	16Fr	8Fr-16Fr	ø0.035"
SPPN0318	18Fr	8Fr-18Fr	ø0.035"
SPPN0320	20Fr	8Fr-20Fr	ø0.035"
SPPN0322	22Fr	8Fr-22Fr	ø0.035"

Set with Drainage Catheter add D (for example SPPN0308+D)

Ureteral Catheters

To inject contrast liquid

With stylet
Made of PEBAX
Available in luer connector



Straight Tip Closed

Product Codes	Length(cm)	Size
SPUC0140	70	4Fr
SPUC0150	70	5Fr
SPUC0160	70	6Fr

Straight Tip Open


Product Codes	Length(cm)	Size
SPUC0240	70	4Fr
SPUC0250	70	5Fr
SPUC0260	70	6Fr

Whistle Tip

Product Codes	Length(cm)	Size
SPUC0340	70	4Fr
SPUC0350	70	5Fr
SPUC0360	70	6Fr

Puncture Needles

Used percutaneous Nephrostomy and
all other percutaneous applications



Product Codes	Size (GA)	Length(cm)
SPND1824	18G	24cm
SPND2024	20G	24cm



Endoscope Protection Cover

Product Codes	Descriptions	Size (Fr)
SPEP0101	Endoscope Protection cover	14×150cm
SPEP0201	Endoscope Protection cover	14×200cm
SPEP0301	Endoscope Protection cover	20×150cm
SPEP0401	Endoscope Protection cover	20×200cm



Surgical Adhesive Membrane

Product Codes	Descriptions	Size (Fr)
SPSF0145	Surgical Adhesive Membrane	45x45cm

Grasping Forceps



Product Codes	Descriptions	Size (Fr)
SPGF0151	Grasping Forceps,rigid	5Fr×600mm
SPGF0251	Grasping Forceps,flexible	5Fr×600mm

Biopsy Forceps



Product Codes	Descriptions	Size (Fr)
SPBF0151	Biopsy forceps,rigid	5Fr×600mm
SPBF0251	Biopsy forceps,flexible	5Fr×600mm

Electrodes

For Prostate Surgery



Sphere Electrode (5mm Ball)



Necklace Electrode



Roller Electrode



Collins Knife



Cutting Loop for 30° Optic



Cold Knife Electrode



Kriloktrode Electrode
(Simultaneous Cutting & Vaporization)



Spike Roller Electrode

Electrodes MONOPOLAR

STORZ, 2 Stem, 30° Degrees, 24Fr

Cutting Loop	SPED1121
Collin's Knife	SPED1122
Roller Electrode	SPED1123
Ball Electrode 3mm	SPED1124

Ball Electrode 5mm	SPED1125
Kriloktrode	SPED1126
Cold Knife	SPED1127
Necklace	SPED1128

Spike 3mm	SPED1129
Spike 5mm	SPED11210

STORZ, 1 Stem, 30° Degrees, 24Fr

Cutting Loop	SPED1111
Collin's Knife	SPED1112
Roller Electrode	SPED1113
Ball Electrode 3mm	SPED1114

Ball Electrode 5mm	SPED1115
Kriloktrode	SPED1116
Cold Knife	SPED1117
Necklace	SPED1118

Spike 3mm	SPED1119
Spike 5mm	SPED11110

OLYMPUS, 1 Stem, 12° Degrees, 24Fr

Cutting Loop	SPED2211
Collin's Knife	SPED2212
Roller Electrode	SPED2213

Ball Electrode 3mm	SPED2214
Sphere Electrode 5mm	SPED2215
Kriloktrode	SPED2216

OLYMPUS, 1 Stem, 30° Degrees, 24Fr

Cutting Loop	SPED2111
Collin's Knife	SPED2112
Roller Electrode	SPED2113

Ball Electrode 3mm	SPED2114
Sphere Electrode 5mm	SPED2115
Kriloktrode	SPED2116

WOLF, 1 Stem, 15° Degrees, 24Fr

Cutting Loop	SPED3311
Collin's Knife	SPED3312
Roller Electrode	SPED3313

Ball Electrode 3mm	SPED3314
Sphere Electrode 5mm	SPED3315
Kriloktrode	SPED3316

Silicone Foley Catheters



2-Way

Product Codes	Balloon Size	Size (Fr)	Length(cm)
SPFC2208	3CC	8Fr	31
SPFC2210	3CC	10Fr	31
SPFC2212	5CC	12Fr	40
SPFC2214	5CC	14Fr	40
SPFC2216	5CC	16Fr	40
SPFC2218	5CC	18Fr	40
SPFC2220	5CC	20Fr	40
SPFC2222	5CC	22Fr	40
SPFC2224	5CC	24Fr	40
SPFC4216	30CC	16Fr	40
SPFC4218	30CC	18Fr	40
SPFC4220	30CC	20Fr	40
SPFC4222	30CC	22Fr	40
SPFC4224	30CC	24Fr	40

3-Way

Product Codes	Balloon Size	Size (Fr)	Length(cm)
SPFC4314	30CC	14Fr	40
SPFC4316	30CC	16Fr	40
SPFC4318	30CC	18Fr	40
SPFC4320	30CC	20Fr	40
SPFC4322	30CC	22Fr	40
SPFC4324	30CC	24Fr	40
SPFC4326	30CC	26Fr	40

Latex Foley Catheters



2-Way

Product Codes	Balloon Size	Size (Fr)	Length(cm)
SPFC1208	3CC	8Fr	27
SPFC1210	3CC	10Fr	27
SPFC1212	5CC	12Fr	40
SPFC1214	5CC	14Fr	40
SPFC1216	5CC	16Fr	40
SPFC1218	5CC	18Fr	40
SPFC1220	5CC	20Fr	40
SPFC1222	5CC	22Fr	40
SPFC1224	5CC	24Fr	40
SPFC1226	5CC	26Fr	40
SPFC1228	5CC	28Fr	40
SPFC1230	5CC	30Fr	40
SPFC3212	30CC	12Fr	40
SPFC3214	30CC	14Fr	40
SPFC3216	30CC	16Fr	40
SPFC3218	30CC	18Fr	40
SPFC3220	30CC	20Fr	40
SPFC3222	30CC	22Fr	40
SPFC3224	30CC	24Fr	40
SPFC3226	30CC	26Fr	40
SPFC3228	30CC	28Fr	40
SPFC3230	30CC	30Fr	40

3-Way

Product Codes	Balloon Size	Size (Fr)	Length(cm)
SPFC3314	30CC	14Fr	40
SPFC3316	30CC	16Fr	40
SPFC3318	30CC	18Fr	40
SPFC3320	30CC	20Fr	40
SPFC3322	30CC	22Fr	40
SPFC3324	30CC	24Fr	40
SPFC3326	30CC	26Fr	40

Urine Meter

400mL capacity meter
with 2500 mL drain bag
60In. of 11/32In. STAR tube
Needless sample port

Product Codes
SPUB0303



Urine Bag

Anti-reflux valve
42 In. of 9/32 In. STAR tube
Needless sample port

Product Codes
SPUB0101



With universal double hanger,
48 In. of 11/32 In. STAR tube,
Anti-reflux drip chamber
Air vent

Product Codes
SPUB0202



Leg Bag

600ml/900ml
Latex free straps
Twist valve and slide-tap
18In. Extension Tubing(optional)



Product Codes	Description	Size (ml)
SPUB0414	T-Tap Valve	600ml
SPUB0415	T-Tap Valve	900ml
SPUB0424	Twist Valve	600ml
SPUB0425	Twist Valve	900ml

Leg Bag with 18 In. Extension Tubing add T (for example SPUB0414+T)

Intermittent Catheter

Polished smooth eyelets
Latex-Free, BPA-Free, DEHP-Free
Straight & Coude Tip
PU & PVC Material

Straight Tip

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC1108	Pediatric	8Fr	10"
SPNC1110	Pediatric	10Fr	10"
SPNC1112	Pediatric	12Fr	10"
SPNC1114	Pediatric	14Fr	10"
SPNC1206	Male	6Fr	16"
SPNC1208	Male	8Fr	16"
SPNC1210	Male	10Fr	16"
SPNC1212	Male	12Fr	16"
SPNC1214	Male	14Fr	16"
SPNC1216	Male	16Fr	16"
SPNC1218	Male	18Fr	16"
SPNC1308	Female	8Fr	7.5"
SPNC1310	Female	10Fr	7.5"
SPNC1312	Female	12Fr	7.5"
SPNC1314	Female	14Fr	7.5"
SPNC1314N	Female	14Fr	6"
SPNC1316	Female	16Fr	7.5"
SPNC1318	Female	18Fr	7.5"

Olive Tip Coude

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC2212	Male	12Fr	16"
SPNC2214	Male	14Fr	16"
SPNC2216	Male	16Fr	16"
SPNC2218	Male	18Fr	16"

Tapered Tip

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC3210	Male	10Fr	16"
SPNC3212	Male	12Fr	16"
SPNC3214	Male	14Fr	16"
SPNC3216	Male	16Fr	16"
SPNC3218	Male	18Fr	16"

Clear & Free Hydrophilic Catheter (Coated Catheter Only)

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC4106	Pediatric	6Fr	10"
SPNC4108	Pediatric	8Fr	10"
SPNC4110	Pediatric	10Fr	10"
SPNC4212	Standard	12Fr	16"
SPNC4214	Standard	14Fr	16"
SPNC4216	Standard	16Fr	16"
SPNC4218	Standard	18Fr	16"
SPNC4310	Female	10Fr	7.5"
SPNC4312	Female	12Fr	7.5"
SPNC4314	Female	14Fr	7.5"
SPNC4316	Female	16Fr	7.5"
SPNC4318	Female	18Fr	7.5"

Clear & Free Hydrophilic Catheter (with Easy Brust Water Pouch only)

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC5106	Pediatric	6Fr	10"
SPNC5108	Pediatric	8Fr	10"
SPNC5110	Pediatric	10Fr	10"
SPNC5212	Standard	12Fr	16"
SPNC5214	Standard	14Fr	16"
SPNC5216	Standard	16Fr	16"
SPNC5218	Standard	18Fr	16"
SPNC5310	Female	10Fr	7.5"
SPNC5312	Female	12Fr	7.5"
SPNC5314	Female	14Fr	7.5"
SPNC5316	Female	16Fr	7.5"
SPNC5318	Female	18Fr	7.5"

Clear & Free Hydrophilic Catheter (with Easy Brust Water Pouch with Touch-Free Sleeve)

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC6108	Pediatric	8Fr	10"
SPNC6110	Pediatric	10Fr	10"
SPNC6212	Standard	12Fr	16"
SPNC6214	Standard	14Fr	16"
SPNC6216	Standard	16Fr	16"
SPNC6218	Standard	18Fr	16"
SPNC6310	Female	10Fr	7.5"
SPNC6312	Female	12Fr	7.5"
SPNC6314	Female	14Fr	7.5"
SPNC6316	Female	16Fr	7.5"
SPNC6318	Female	18Fr	7.5"

Clear & Free Hydrophilic Catheter-Coude (with Easy Brust Water Pouch and Touch-Free Sleeve)

Product Codes	Descriptions	Size (Fr)	Length(In)
SPNC7210	Coude-Standard	10Fr	16"
SPNC7212	Coude-Standard	12Fr	16"
SPNC7214	Coude-Standard	14Fr	16"
SPNC7216	Coude-Standard	16Fr	16"