

ESSx™ Microdebrider

**Bringing power and efficiency to
your operating room.**



**Introducing the new,
ergonomically advanced ESSx shaver.**

ESSxTM Microdebrider

stryker[®]


Neuro Spine ENT

ESSx Microdebrider

5290-601-200



- Max RPM in oscillation mode (standard): 5,000 RPM
- Improved ergonomics and balance
- Increased cutting efficiency
- Straight through suction: Helps ensure less clogging
- Improved patented flush feature: Surgeon can manually unlog hand piece using its own irrigation out flow allowing for quicker removal of clogs and decreased OR time
- Integrated scope cleaning: Cleans the endoscope allowing for less manual cleaning during the case



ENT Cutters and Burs

- All blades are irrigated through the entire channel including the tip
- Preserves mucous membranes



SharpEdge

Dia.	Order #	Length
4.0 mm	5290-340-000	8cm

Aggressive


Dia.	Order #	Length
2.5 mm	5290-628-000	8cm
3.5 mm	5290-638-000	8cm
4.0 mm	5290-648-000	8cm
4.0 mm	5290-011-000	11cm

Aggressive Serrated

Dia.	Order #	Length
3.5 mm	5290-635-000	8cm
4.0 mm	5290-645-000	8cm

Jaguar

Dia.	Order #	Length
3.5 mm	5290-530-000	8cm
4.0 mm	5290-540-000	8cm



Angled Aggressive 40°

Dia.	Order #	Length
4.0 mm	5290-740-100	11cm

Angled Aggressive 40° Convex

Dia.	Order #	Length
4.0 mm	5290-740-200	11cm

Angled Aggressive 12°

Dia.	Order #	Length
4.0 mm	5290-744-100	11cm

Hooded Round Bur

Dia.	Order #	Length
3.0 mm	5290-647-000	8cm

Unhooded Round Bur

Dia.	Order #	Length
2.0 mm	5290-220-000	8cm

Accessories



iSwitch Receiver
277-200-100

CORE Console
5400-050-000

Irrigation Cassette
5400-050-001



iSwitch Footswitch
277-100-100

Bi-directional Foot Pedal
5100-008-000



Mini-Trephination Set
0290-222-000

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker. Products referenced with [®] designation are registered trademarks of Stryker.

Literature Number: 9100-001-091 Rev. None
UnDe/PS

Copyright © 2006 Stryker
Printed in USA

Stryker[®] Neuro Spine ENT
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210

www.stryker.com