




Always-On Tip Tracked[®] Instrument Catalog


SPiN THORACIC NAVIGATION SYSTEM™




Superior Navigation

Our revolutionary Always-On Tip Tracked® Instruments feature electromagnetic sensors that enable you to see straight to the target throughout the entire procedure—even during biopsy, without the need for fluoroscopy.

INS-5900
SPiN ACCESS CATHETER™
KIT (90 degree catheter)



Also available in 180 degree catheter
kit (INS-5905)



[+] SPiN Access Catheter™ fits down therapeutic scope with 2.8mm WC or larger



Always-On Tip Tracked® Pulmonary Products

PART NUMBER	PRODUCT INFORMATION
INS-0048	vPad™ Connector Cable
INS-0050	vPad™ Patient Tracker (Box of 5)
INS-0304	Always-On Tip Tracked® Guidewire, 1.0mm, 735mm L (Box of 5)
INS-0305	Always-On Tip Tracked® Guidewire, 1.0mm, 950mm L (Box of 5)
INS-0322	View Peripheral Catheter™, 3.2mm OD, 2.0mm WC (Box of 5)
INS-5440	View Peripheral Catheter™, 3.8mm OD, 2.5mm WC (Box of 5)
INS-0352	Always-On Tip Tracked® Brush, 15mm L, 1.8mm OD (Box of 5)
INS-0362	Always-On Tip Tracked® Forceps, 1.8mm OD - Oval Cup (Box of 5)
INS-0372	Always-On Tip Tracked® Forceps, 1.8mm OD - Serrated Cup (Box of 5)
INS-0382	Always-On Tip Tracked® 19ga ANSO Histology Needle, 1.8mm OD (Box of 5)
INS-0392	Always-On Tip Tracked® 21ga ANSO Cytology Needle, 1.8mm OD(Box of 5)
INS-5410	Always-On Tip Tracked® 22ga SPiN Flex™ ANSO Nitinol Needle, 12mm L, 1.8mm OD (Box of 5)
INS-5411	0.016" Stylet, 1235mm L (Box of 10)
INS-5500	View Optical Probe* (Box of 5)
INS-5600	SPiN Perc® Kit (Box of 5)
INS-5620	SPiN Perc® Kit, 20cm (Box of 5)
INS-5700	SPiN EWC™ Biopsy Guide Kit (Box of 5)
INS-5900	SPiN Access Catheter™ 90 Degree for use with Olympus Scopes (Box of 5)
INS-5905	SPiN Access Catheter™ 180 Degree for use with Olympus Scopes (Box of 5)
INS-5450	SPiN Xtend™ (Box of 5)
INS-5300	Always-On Tip Tracked® Triple Needle Brush, 12 mm L, 1.8mm OD (Box of 5)
INS-5610	SPiN Perc® Localization Needle, 1cm (Box of 5)
INS-5611	SPiN Perc® Localization Needle, 2cm (Box of 5)

INS-5600
SPiN PERC® KIT


19ga x 105mm Always-On Tip Tracked® biopsy needle,
20ga x 15cm biopsy gun, and 20ga x 15cm FNA needle.

Also available in 20cm Kit: 19ga x 155mm Always-On Tip
Tracked® biopsy needle, 20ga x 20cm biopsy gun, and
20ga x 20cm FNA needle (INS-5620).





Extend reach to targets up to
4cm outside airways



**INS-5450
SPIN XTEND™**

Always-On Tip Tracked®
4cm Extendable Needle
21ga Nitinol Needle




**INS-5410
SPIN FLEX™ NEEDLE**
22ga, 12mm L, 1.8mm OD

Laser-etched
Nitinol Needle



Steer 0° to 180°
in all directions



**INS-5440
VIEW PERIPHERAL CATHETER™**
3.8mm OD, 2.5mm WC

Also available in 3.2mm OD, 2.0mm WC
(INS-0322)

**INS-0382
ALWAYS-ON TIP TRACKED® NEEDLE**

19ga, 12mm L, 1.8mm OD
Also available in 21ga, 12mm L, 1.8mm OD
(INS-0392)



**INS-0352
ALWAYS-ON TIP TRACKED®
BRUSH**

15mm L, 1.8mm OD




**INS-5300
ALWAYS-ON TIP TRACKED®
TRIPLE NEEDLE BRUSH**

12mm L, 1.8mm OD



**INS-0304
ALWAYS-ON TIP
TRACKED® GUIDEWIRE**

1.0mm, 735mm L
Also available in 1.0mm,
950mm L (INS-0305)



**INS-5500
VIEW OPTICAL
PROBE***




The View Optical
Probe fits down
working channel of
the View Peripheral
Catheter™




**INS-0362
ALWAYS-ON TIP
TRACKED® FORCEPS**

Oval Cup, 1.8mm OD

Also available in Serrated
Cup, 1.8mm OD (INS-0372)



INS-0050
vPad™
Patient Tracker




INS-0048
vPad™
Connector Cable




Always-On Tip Tracked® Radiology Products

PART NUMBER	PRODUCT INFORMATION
INS-5017	17ga x 105mm Tip Tracked Bx Needle (Box of 5)
INS-5024	17ga x 150mm Tip Tracked Bx Needle (Box of 5)
INS-5028	19ga x 155mm UTW Tip Tracked Bx Needle, Chiba Point (Box of 5)
INS-5029	19ga x 105mm UTW Tip Tracked Bx Needle, Chiba Point (Box of 5)
INS-5032	19ga x 150mm UTW Tip Tracked Bx Needle (Box of 5)
INS-5033	19ga x 60mm UTW Tip Tracked Bx Needle (Box of 5)
INS-5034	17ga x 123mm Tip Tracked Bx Needle Lg, Chiba Point (Box of 5)
INS-5036	22ga x 125mm UTW Tip Tracked Bx Needle (Box of 5)
INS-5039	Always-On vTrack Universal Tracker, SM/CL, 12-18ga (Box of 5)
INS-5040	Always-On vTrack Universal Tracker, LG 7-12ga (Box of 5)

INS-5610
SPiN PERC® LOCALIZATION NEEDLE
19ga needle with 1cm marking zone



INS-5611
SPiN PERC® LOCALIZATION NEEDLE
19ga needle with 2cm marking zone



**TIP TRACKED
BIOPSY NEEDLES**
Available in 17ga, 19ga &
22ga configurations and
multiple lengths



INS-5039
**ALWAYS-ON vTRACK
UNIVERSAL TRACKER**
SM/CL, 12-18ga



For tracking rigid
ablation probes



INS-5040
**ALWAYS-ON vTRACK
UNIVERSAL TRACKER**
LG, 7-12ga



For tracking rigid
ablation probes



Need help finding or ordering an instrument? Contact Us

Tell us how we can help by sending an email to orders@veranmedical.com

Veran Medical Technologies

1908 Innerbelt Business Center Dr. | St. Louis, MO 63114

P 314.659.8500 | **F** 314.659.8560 | **E** orders@veranmedical.com

* Used with the SPIN View System. CAUTION - CLASS 4 LASER RADIATION WHEN OPEN AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION. VISIBLE LASER RADIATION AVOID DIRECT EYE EXPOSURE CLASS 3R LASER PRODUCT

