




PRODUCTS FOR PERCUTANEOUS TRACHEOSTOMY

Arndt Multiport Airway Adapter

Permits simultaneous bronchoscopy and ventilation of the airway.



The Arndt Multiport Adapter and Cook Airway Exchange Catheter can be used during balloon-assisted tracheostomy to provide:

- **Continuous visualization** to facilitate precise placement of the introducer and help prevent posterior tracheal wall damage.
- **Continuous airway access** via the Cook Airway Exchange Catheter, which serves as a guide for reintubation and can deliver oxygen during the procedure if necessary.

Global
Product
Number

Order
Number

Set

G44105

C-CAE-11.0-100-AMPA

Weinmann Tracheostomy Exchange Set

Facilitates smooth and efficient redilatation of the stoma if necessary for tracheostomy tube exchange.



Global Product Number	Order Number	Loading Dilator Fr
Set		
G31447	C-PTIS-100-HC-EXCH	26.0, 28.0



www.cookmedical.com

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 USA
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010 CC-BM-PTBSA-PDT-EN-201010