



Head Lights ●

Fibreoptic Cold Light Head Lights		REF
Cold light headlamp (standard) Complete with leather-upholstered headband, 2.2 m fibreoptic light cable installed on right, luminous field 50 to 210 mm 0° (straight) cable head	006-08315	
Cold light headlamps (standard) Complete with leather-upholstered headband, 2.2 m fibreoptic light cable installed on right, luminous field 50 to 210 mm 90° (angled) cable head	006-08317	
Microspot surgical cold light headlamp Focusable fibreoptic headlamp with high-performance optics and a light, comfortable headband. For the illumination of the smallest bodily cavities. Illuminous field of 7 to 63 mm. 0° (straight) cable head	006-08370	
Microspot surgical cold light headlamp Focusable fibreoptic headlamp with high-performance optics and a light, comfortable headband. For the illumination of the smallest bodily cavities. Illuminous field of 7 to 63 mm. 90° (angled) cable head	006-08372	
Warm Light Head Light		
CLAR-S headlamps Headlamp complete with a Ø 55 mm concave mirror, adjustable, ergonomic headband and 2 m cable with banana plugs	006-0812106	

LED Cold Light Head Light	REF	
<p>Standard surgical LED headlamp (black with headband)</p> <p>High-power, focusable headlamp with a white, high-performance LED. The Li-Ion battery integrated in the headband makes this light, comfortable model especially suitable for mobile applications.</p> <p>Total weight: 350 g (with battery pack on the headband) Charging time: Approx. 4 hours from completely dead Operating period: Approx. 3.5 hours Working distance: 20 to 60 cm Illuminous field diameter: 11 to 130 mm Colour temperature: Approx. 6,000° K</p>	006-08795-S	
<p>Standard surgical LED headlamp (white with headband)</p> <p>High-power, focusable headlamp with a white, high-performance LED. The Li-Ion battery integrated in the headband makes this light, comfortable model especially suitable for mobile applications.</p> <p>Total weight: 350 g (with battery pack on the headband) Charging time: Approx. 4 hours from completely dead Operating period: Approx. 3.5 hours Working distance: 20 to 60 cm Illuminous field diameter: 11 to 130 mm Colour temperature: Approx. 6,000° K</p>	006-08795-W	
<p>Profi 3S LED headlamp</p> <p>Extra bright LED headlamp with adjustable light intensity: Weight: 336 g, operating period: 10 hours (mPack), brightness: 50,000 lux, colour temperature: 6,250° K, luminous field size: 30-80 mm (at working distance of 420 mm)</p> <p>With mPack (battery) and plug-in transformer</p> <p>Binocular loupes for the 3S LED headlamp available on request</p>	122-00831313	