

Instrument Cabinets

Are you looking for stable, spacious and mobile instrument cabinets or charts? We offer a wide range of instrument cabinets and video charts, which not only have the same design as the MODULA treatment units, but which you can also customize to suit your personal requirements.

All the cabinets and trolleys are composed of a stable, electrolytically galvanised plate housing with powder-coated stoved enamelling. The cabinets are all mobile and equipped with an acrylic glass cover. We offer you every model with basic equipment that you can expand to suit your requirements. Not only can you select the individual components in and on the cabinet however you wish, but you can also even choose the colour of the cabinet (20 RAL colours).

HANSEAT-L

REF: 0300-10000

- High-quality, electrolytically galvanised sheet steel housing
- Powder-coated stoved enamelling in RAL colours (see standard colour table)
- Mobile unit with stainless steel base, antistatic heavy-duty castors
- 4 fully extendable and soft closing drawers
- 2 instrument levels with stainless steel trays
- **Melamine instrument tray for used instruments**
- **Waste bin opened with foot pedal, closes quietly automatically**
- Acrylic glass instrument cover, sinkable behind the cabinet
- Pull-out writing/preparation desk with stainless steel surface



HANSEAT-XL

REF: 0300-10100

- High-quality, electrolytically galvanised sheet steel housing
- Powder-coated stoved enamelling in RAL colours (see standard colour table)
- Mobile unit with stainless steel base, antistatic heavy-duty castors
- 4 continuous drawers, fully extending with soft close system
- 2 instrument levels with stainless steel trays
- Acrylic glass instrument cover, sinkable behind the cabinet
- Pull-out writing/preparation desk with stainless steel surface



HANSEAT-M

REF: 0300-20000

- High-quality, electrolytically galvanised sheet steel housing
- Powder-coated stoved enamelling in RAL colours (see standard colour table)
- Mobile unit with stainless steel base, antistatic heavy-duty castors
- 4 fully extendable and soft closing drawers
- 2 instrument levels with stainless steel trays
- Acrylic glass instrument cover, sinkable behind the cabinet



| Options for MODULA Instrument Cabinets | REF | HANSEAT-L | HANSEAT-XL | HANSEAT-M |
|---|------------|-----------|------------|-----------|
| Housing Options | | | | |
| LED instrument illumination with changing RGB colours (red, green, blue) incl. remote control | 0110-20010 | ● | ● | ● |
| 3 rd instrument level with extendable central level (not together with REF: 0110-20025) | 0110-20000 | ● | ● | ● |
| Heatable instrument storage surface (only for the top level respectively) | 0110-20020 | ● | ● | ● |
| Mirror warming block for approx. 50 mirrors (not together with REF: 0110-20025 or 0110-20025) | 0110-20110 | ● | ● | ● |
| Storage compartment for medicament bottles, integrated in the instrument level (not together with REF: 0110-20110) | 0110-20025 | ● | ● | ● |
| Upper shelf above instrument level with powder-coated aluminium shelf (integrated board surface for devices) | 0110-20060 | ● | ● | ● |
| Upper shelf above instrument level with matt acrylic glass shelf (integrated board surface for devices) | 0110-20065 | ● | ● | ● |
| Lower shelf (instead of drawer) | 0110-21000 | ● | - | ● |
| Flat, half-open drawer for instruments | 0110-21300 | ● | - | ● |
| Storage compartment with door (instead of waste bin) | 0110-21200 | ● | - | - |
| Cotton wool and tongue swab dispenser (attached on right side of unit) | 0110-20130 | ● | ● | ● |
| Cotton wool and tongue swab dispenser (integrated in the unit) | 0110-20135 | ● | - | ● |
| Desk extension S (approx. 50 cm extension) | 0110-20150 | ● | ● | - |
| Desk extension M (approx. 100 cm extension) | 0110-20160 | ● | ● | ● |
| Desk extension M with sliding keyboard shelf (approx. 100 cm extension) | 0110-20165 | ● | ● | ● |
| Desk extension L (approx. 120 cm extension) | 0110-20170 | ● | ● | ● |
| Desk extension L with sliding keyboard shelf (approx. 120 cm extension) | 0110-20175 | ● | ● | ● |
| NEW: UV-C Hygiene Drawer (instead of standard drawer) | 0110-91000 | ● | - | - |

| Options for MODULA Instrument Cabinets | REF | HANSEAT-L | HANSEAT-XL | HANSEAT-M |
|--|------------|-----------|------------|-----------|
| Accessories Adaptation | | | | |
| Holder for camera head (magnetic) | 0110-20082 | ● | ● | ● |
| Stand column for swivel arms (installation on treatment unit) | 0110-71030 | ● | ● | - |
| Stand column + monitor holder for 10" to 19" TFT monitor (installation on unit)* | 0110-71035 | ● | ● | - |
| Swivel arm with storage board for additional devices (Installed on microscope column or stand column)* | 0110-71010 | ● | ● | - |
| Swivel arm with storage board for endoscopy camera "starCAM" or "swift-CAM" (installation on microscope column or stand column)* | 0110-71041 | ● | ● | - |
| Swivel arm with storage board for endoscopy camera "XE-50/TFT" (installation on microscope column or stand column)* | 0110-71042 | ● | ● | - |
| Swivel arm with storage board for diagnostic devices (installation on microscope column or stand column)* | 0110-71043 | ● | ● | - |
| Instrument Storage | | | | |
| Stainless steel instrument tray attachment E1 for top level (for units with REF: 0110-20025) | 0160-10040 | - | - | ● |
| Stainless steel instrument tray attachment E2 for top level | 0160-10050 | - | - | ● |
| Stainless steel instrument tray attachment E3 for lower/middle level | 0160-10060 | - | - | ● |
| Stainless steel instrument tray attachment E4 for top level (for units with REF: 0110-20000) | 0160-10070 | - | - | ● |
| Stainless steel instrument tray attachment E5 for top level (for units with REF: 0110-20025) | 0160-10000 | ● | ● | - |
| Stainless steel instrument tray attachment E6 for top level | 0160-10010 | ● | ● | - |
| Stainless steel instrument tray attachment E7 for lower/middle level | 0160-10020 | ● | ● | - |
| Stainless steel instrument tray attachment E8 for top level (for units with REF: 0110-20000) | 0160-10030 | ● | ● | - |
| Stainless steel ear speculum holder, small (3x6) | 0160-21000 | ● | ● | ● |
| Stainless steel ear speculum holder, medium (4x6) | 0160-21010 | ● | ● | ● |
| Stainless steel tongue depressor holder, small (3 rows) | 0160-21200 | ● | ● | ● |
| Stainless steel tongue depressor holder, medium (4 rows) | 0160-21210 | ● | ● | ● |

*Only one swivel arm per stand or microscope column possible

| Options for MODULA Instrument Cabinets | REF | HANSEAT-L | HANSEAT-XL | HANSEAT-M |
|--|------------|-----------|------------|-----------|
| Instrument Storage | | | | |
| Stainless steel instrument tray for used instruments | 0160-12000 | - | - | ● |
| Silicone insert for instrument tray | 0160-11020 | - | - | ● |
| Stainless steel cover for instrument tray | 0160-12100 | - | - | ● |
| Splash protection cover for instrument tray | 0160-12400 | - | - | ● |
| Melamine instrument tray for used instruments | 0160-11000 | ● | ● | - |
| Silicone insert for instrument tray | 0160-11030 | ● | ● | - |
| Stainless steel basket for instrument tray | 0160-11010 | ● | ● | - |
| Splash protection cover for instrument tray | 0160-12300 | ● | ● | - |
| Stainless steel spray bottle holder with 3 drug nebulisers (integration in the bottom level) | 0160-21300 | ● | ● | ● |

Device/video trolleys:

- High-quality, electrolytically galvanised sheet steel housing
- Powder-coated stoved enamelling in RAL colours (see standard colour table)
- Mobile unit on stainless steel base
- 1 drawer
- Double-walled storage shelves, 480 x 480 mm
- The storage shelves can be adjusted in increments of 50 mm
- The intermediate shelves are very strong with a load-bearing capacity of approx. 50 kg
- Carrier sections each with 4 openings 50 mm in diameter
- Perforated sheet rear panel with 4 cable holes 50 mm in diameter
- 4 double-caster rubber wheels, the front two with lockable brake



| Ordering Information | REF |
|---|------------|
| Device/video trolley L 900, 4 levels (WxHxD: 900 x 555 x 450 mm) | 0300-70000 |
| Device/video trolley H 1.20 4 levels (WxHxD: 1,200 x 555 x 450 mm) | 0300-60000 |

Additional equipment available on request