

LIGHT SOURCES

L9000 Light Source



L9000 LED Light Source
LED technology

- » No bulb replacement
- » Designed to produce less environmental waste
- » Safelight cable compatible

- 220-210-000** L9000 LED Light Source
- 6000-001-020** Firewire Cable
Not included with L9000, but needed for light source control from the 1288 camera head

X8000 Light Source



- 220-200-000** X8000 Xenon Light Source

Light Source Accessories

Safelight Fiberoptic Light Cables and Adapters

A coiled white fiberoptic light cable with a green connector on one end and a blue connector on the other.

- 233-050-199** Safelight Kit:
Safelight Cable (233-050-100) & Safelight Scope Adapter (233-050-087)
- 233-050-100** 5.0 mm x 10 ft Clear-Case Safelight Fiberoptic Light Cable
- 233-050-200** 5.0 mm x 10 ft Opaque-Case Safelight Fiberoptic Light Cable
- 233-050-087** Safelight Henke/Dyonics/Storz Scope Adapter
- 233-050-088** Safelight Wolf Scope Adapter
- 233-050-086** Safelight Olympus Scope Adapter
- 233-050-192** Safelight Circon Scope Adapter



233-065-010
6.5 mm x 10 ft



233-050-090
5.0 mm x 10 ft



233-050-064
5.0 mm x 10 ft
Clear Case



233-050-065
5.0 mm x 10 ft



233-050-069
5.0 mm x 10 ft
Clear Case



233-050-067
2.0 mm x 7.5 ft
Clear Case




233-050-084
5.0 mm x 10 ft
Non-Transparent

Fiberoptic Light Cables

- 233-065-010** 6.5 mm x 10 ft Fiberoptic Light Cable
Wolf Scope & Stryker Light Source Fittings
- 233-050-090** 5.0 mm x 10 ft Fiberoptic Light Cable
Universal Scope and Stryker Light Source Fittings
- 233-050-064** 5.0 mm x 10 ft Clear Case Fiberoptic Light Cable
Storz Scope & Stryker Light Source Fittings, Screw-On
- 233-050-065** 5.0 mm x 10 ft Fiberoptic Light Cable
Stryker Scope & Stryker Light Source Fittings
- 233-050-069** 5.0 mm x 10 ft Clear Case Fiberoptic Light Cable
Wolf Scope & Stryker Light Source Fittings, Snap-On
- 233-050-067** 2.0 mm x 7.5 ft Clear Case Fiberoptic Light Cable
Stryker Scope & Stryker Light Source Fittings
- 233-050-084** 5.0 mm x 10 ft Non-Transparent Fiberoptic Cable
Storz Scope and Stryker Light Source Fittings, Screw-On

Replacement Bulbs



- 220-201-000** Replacement Bulb for X8000 (300 Watts)
- 220-191-000** Replacement Bulb for X7000 (300 Watts)

Light Source Adapters

- » Connects to Light Source
- » For use with universal fiberoptic light cable, 233-050-090

- 233-050-071** Stryker/ACMI/Zimmer Light Source Adapter
- 233-050-072** Olympus Light Source Adapter
- 233-050-073** Storz Light Source Adapter
- 233-050-074** Wolf/Dyonics Light Source Adapter



233-050-075



233-050-076



233-050-077



233-060-000

Scope Adapters

- » Connect to Scope
- » For use with universal fiberoptic light cable, 233-050-090

233-050-075	Stryker/ACMI Scope Adapter
233-050-076	Wolf/Dyonics Scope Adapter
233-050-077	Storz/Olympus Scope Adapter
233-060-000	ACMI Universal Scope Adapter



233-080-000



233-081-000



233-082-000

Scope Adapters for Non-Autoclavable Scopes

- » Black Eyepiece
- » Connect to Scope

233-080-000	Stryker Adapter
233-081-000	Wolf/Dyonics Adapter
233-082-000	Storz Adapter



331-000-100



331-000-200



331-000-300



331-000-400

Scope Adapters for Autoclavable Scopes

- (Blue & Red Eyepiece)
- Connect to Scope

331-000-100	Stryker Adapter
331-000-200	Wolf/Dyonics Adapter
331-000-300	Storz Adapter
331-000-400	Right Angle Adapter