

ENT Microscopes



SOM 22

SOM 62

SOM model range

Standard equipment:

- 3-level magnification changer
- Lens $f = 200$ mm
- Binocular straight tube $f = 125$ mm
- Eyepieces for spectacle wearers
- Illumination unit 15 V, 150 VA

● Examples of additional equipment



**Microscope holder
with fine focussing**



**Beam splitter with
endoscope adapter**



Lenses



Rapid lens changer



Handles

Ordering Information	REF
KAPS SOM 22 (installation in unit) with cold light illumination	535-0500059
KAPS SOM 22 (installation in unit) with LED cold light illumination	
KAPS SOM 32 (wall installation) with cold light illumination	535-0500011
KAPS SOM 62 (floor installation) with cold light illumination	535-0500014
Beam splitter with endoscope adapter	535-101890000
Lens f = 400 mm	535-101670000
Vario lens f = 200-400 mm	535-162700000
Rapid lens changer for 3x changer (without lens)	535-102181000

Additional equipment available on request