



## OPMI pico

Comprising the OPMI pico with 5x changer, integrated green and orange filter, a T-handle and integrated cold light fibreoptic illumination with rapid lamp changer

- Straight wide-field tube 45°  $f = 170$  mm
- 2 eyepieces with 10x magnification
- Lens  $f = 250$  mm
- Wall, ceiling or floor stand including illumination unit and fibreoptic light cable
- 2 sleeves for tubes with screw thread
- Dust cover
- 6-pack of asepsis caps 12 mm





Ordering Information	REF
Zeiss Opmi pico (installation on unit)	547-0500074
Endoport (beam splitter with endoscope adapter)	547-1153971