


Kröber O2

The quietest Oxygen Concentrator

- modern design
- 30.000 operating hours warranty



- simple operation



- easy to maintain



- innovative Details



- coarse dust filter

Kröber O2

Sauerstoffkonzentrator

Oxygen Concentrator

- modular components set up
- easily accessible filter and fuses
- USB – interface
- 6L/min flow
- digitalised flow adjustment
- modern design and simple operation
- new, innovative measuring sensor
- integrated flow adjustment
- 30.000 hours warranty
- whisper silent
- microprocessor-controlled
- favourable price
- Made in Germany

Accessories

- instructions for use
- angular connector for humidifier
- humidifier, refillable
- nasal cannula 2 m and 5 m



Technical data

MDD classification	II a
Operating power	230 V 50 Hz
Operating temperature	operation +5° to +40°C storage -25° to +70°C
Sound level	35 dB (A)
Power consumption	350 W
Dust filter	in rear of device
Fine filter	behind the maintenance lid
Fuses	mains 230 V 2xTT2,5 A H 250 V internal 1xT1,0 A L 250 V
Alarms	temperature, power failure, flow, oxygen, system
Weight	19,8 kg
Dimensions (HxLxW)	(535x203x520) mm without rollers
Manufacturer's warranty	30.000 operating hours (max. 5 years) according to terms of warranty
Flow adjustment	0 to 2 l/min in 0,1 steps 2 to 4 l/min in 0,2 steps 4 to 6 l/min in 0,5 steps
O2-concentration	1 up to 4 l/min 95% -3% 4 up to 5 l/min 90% +/-3% 5 up to 6 l/min 85% +/-3%



www.easmed.com
customer@easmed.com

Singapore
Tel : (65) 6336 0806
Fax : (65) 6338 4132

Sales Representative
Jessica - (65) 9662 3908
Lyn - (65) 8738 7730

Kröber Medizintechnik GmbH

Salzheck 4 • D-56332 Dieblich

Tel. +49 (0) 26 07 / 9 40 40 • Fax +49 (0) 26 07 / 9 40 422

E-Mail info@kroeber.de • Internet www.kroeber.de