

E.N.T.

RIGID AND FLEXIBLE  
ENDOSCOPES

ENDOSCOPES

SINUSCOPES / LARYNGOSCOPES / OTOSCOPE



We offer exclusive design,  
excellent quality and cost  
effective endoscopes.  
Ecleris innovation value

## ENDOSCOPES

### SINUSCOPES / LARYNGOSCOPES / OTOSCOPE

Outstanding image quality.  
Exceptional luminosity and brightness.  
Storz / Wolf and ACMI connections.



SINUSCOPES			
DIAMETER	3.0mm		
WORKING LENGTH	Approx. 140mm		
LOCKING DEVICE	Storz Type		
DIRECTION OF VIEW	0°	30°	70°



SINUSCOPES			
DIAMETER	4.0mm		
WORKING LENGTH	Approx. 175mm		
LOCKING DEVICE	Storz Type		
DIRECTION OF VIEW	0°	30°	45° 70°



OTOSCOPE	
DIAMETER	3.0mm
WORKING LENGTH	Approx. 100mm
LOCKING DEVICE	Storz Type
DIRECTION OF VIEW	0°



LARYNGOSCOPES		
DIAMETER	6.0mm	8.0mm
WORKING LENGTH	Approx. 170mm	Approx. 170mm
LOCKING DEVICE	-	-
DIRECTION OF VIEW	70°	70°



FLEXIBLE SCOPE	
FIELD OF VIEW	85°
DEPTH OF FIELD	2,5 - 50 mm
OUTER DIAMETER INSERTION TUBE	3,6 mm
OUTER DIAMETER DISTAL END	3.4 mm
WORKING LENGTH	300 mm
BENDING RANGE OF DISTAL TIP	up 130° down 130°

Rev. 03.01.13



[www.ecleris.com](http://www.ecleris.com)

**ECLERIS**   
INNOVATION VALUE