


Discover comprehensive products for minimally invasive drainage procedures.



PRODUCTS FOR CENTESIS AND DRAINAGE



Contents

Pneumothorax catheters


Wayne Pneumothorax Catheter - Seldinger Set.....	4
Wayne Pneumothorax Catheter - Trocar Set.....	5
Cook Emergency Pneumothorax Set	6
Pneumothorax Set.....	7
Catheter Aspiration Set for Simple Pneumothorax.....	8
Simple Pneumothorax Aspiration Accessory Set	9

Multipurpose catheters

Thal-Quick Chest Tube Set	10
Thal-Quick Chest Tube Adapter.....	11
Fuhrman Pleural/Pneumopericardial Drainage Set	12
Wood Pleural/Pneumopericardial Drainage Set.....	13
Thoracentesis Set	14

Pericardiocentesis catheters

Pericardiocentesis Set-Straight Catheter.....	15
Pericardiocentesis Set-Pigtail Catheter.....	16-17
Lock Pericardiocentesis Set	18
Teitel Pericardiocentesis Set.....	19



Wayne Pneumothorax Catheter Set - Seldinger Set

Used for the relief of simple, spontaneous, iatrogenic, and tension pneumothorax.



Features and benefits

- The catheter can be used for either Seldinger or trocar placement.
- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The enlarged sideports are designed to enhance the drainage of air.
- The radiopaque catheter material enhances x-ray visualisation.
- Pigtail catheters have a comparable efficacy to chest tubes in patients with pneumothorax.¹


1.Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma Acute Care Surg.* 2011;71(5):1104-1107.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire gage	Wire Guide Length cm	Dilator Fr
Seldinger set									
G56535	C-UTPT-1400-WAYNE-112497-IMH	14	29	19	18	5,9	.038	70	14


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




ECHOTIP® ECHOGENIC ONE-PART PERCUTANEOUS ENTRY NEEDLE




ONE-PART PERCUTANEOUS ENTRY NEEDLE




DILATOR




#11 THUMB SCALPEL



COOK CHEST DRAIN VALVE



NEEDLE HOLDER CUP






VINYL CONNECTING TUBE


Wayne Pneumothorax Catheter Set - Trocar Set

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm
Trocar set						
G56532	C-UTPT-1020-WAYNE-IMH	10.2	29	20	19	16.5
G56533	C-UTPT-1400-WAYNE-IMH	14	29	20	19	16.5


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




#11 THUMB SCALPEL



COOK CHEST DRAIN VALVE



NEEDLE HOLDER CUP



VINYL CONNECTING TUBE



Cook Emergency Pneumothorax Set


Used for emergency relief and temporary management of suspected tension pneumothorax.

Features and benefits

- The reinforced catheter is designed to prevent collapse after placement.
- The set comes with the Cook Chest Drain Valve, which is designed to permit continued re-expansion of the lung, preventing inadvertent backflow of air into the pleural space.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Needle gage
G27292	C-TPTSJ-8.5-6.0-FSNS	8.5	6	15
G55525	C-TPTSJ-8.5-8.7-FSNS	8.5	8.7	18


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




EMERGENCY PNEUMOTHORAX CATHETER NEEDLE ASSEMBLY



12 mL MONOJECT® SYRINGE



COOK CHEST DRAIN VALVE




CATHETER RETENTION DISC



TY-RAP CABLE TIE



VINYL CONNECTING TUBE



Pneumothorax Set


Used to treat spontaneous and tension pneumothorax.

Features and benefits


- The catheter is used for trocar placement.
- The set comes with the Cook Chest Drain Valve, which is designed to permit continued re-expansion of the lung, preventing inadvertent backflow of air into the pleural space
- The radiopaque catheter material enhances x-ray visualisation.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm
Peters set						
G05299	C-TPT-100-PTR	9	20	20	18	20
Paediatric set						
G03465	C-TPT-100	6.3	19	6	21	12
Standard set						
G03301	C-TPT-100	9	29	20	18	20
G03302	C-TPT-200	10.2	29	20	18	20


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




PNEUMOTHORAX CATHETER/NEEDLE ASSEMBLY



COOK CHEST DRAIN VALVE
(PAEDIATRIC AND STANDARD SETS ONLY)



VINYL CONNECTING TUBE



Catheter Aspiration Set for Simple Pneumothorax


Used for the relief of simple, spontaneous or iatrogenic pneumothorax and tension pneumothorax.

Features and benefits


- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The radiopaque catheter material enhances x-ray visualisation.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Needle gage	Needle Length cm
G05296	C-CASP-901-FORD	9	25	10	.038	50	18	7

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




ASPIRATION CATHETER



ONE-PART PERCUTANEOUS ENTRY NEEDLE



PTFE-COATED FIXED CORE STRAIGHT WIRE GUIDE



THREE-WAY CLEAR PLASTIC STOPCOCK



Simple Pneumothorax Aspiration Accessory Set


Used for the relief of simple, spontaneous or iatrogenic pneumothorax and tension pneumothorax.

Features and benefits


- The Cook Chest Drain Valve is designed to permit continued re-expansion of the lung, preventing inadvertent backflow of air into the pleural space.
- The Cook Chest Drain Valve can be used for the transportation of patients who need to have pleural drainage catheters in place for long durations, or for outpatient procedures.

Order Number	Reference Part Number	Connecting Tube Fr	Connecting Tube Length cm
G05297	C-CASP-A-FORD	14	30

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.



COOK CHEST DRAIN VALVE



VINYL CONNECTING TUBE



Thal-Quick Chest Tube Set


Used for the percutaneous introduction of a chest tube for pleural fluid drainage.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The centimetre-marked dilators allow for additional control over the insertion of a chest tube into the pleural space.
- The sideports are positioned at the distal end of the chest tube.
- The radiopaque catheter material enhances x-ray visualisation.
- The double lumen set has an additional lumen available for aspiration or infusion.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
Single lumen set									
G06786	C-TQTS-800	8	18	3	18	5	.038	55	9
G06787	C-TQTS-1000	10	20	3	18	5	.038	55	9, 11
G05161	C-TQTS-1200	12	22	3	18	7	.038	55	9, 13
G05777	C-TQTS-1400	14	41	4	18	7	.038	65	9, 16
G04891	C-TQTS-1600	16	41	4	18	7	.038	65	9, 13, 18
G05453	C-TQTS-1800	18	41	4	18	7	.038	65	9, 14, 20
G05454	C-TQTS-2000	20	41	4	18	7	.038	65	9, 14, 20
G04137	C-TQTS-2400	24	41	4	18	7	.038	65	9, 18, 26
G04139	C-TQTS-2800	28	41	4	18	7	.038	65	11, 20, 30
G07243	C-TQTS-3200	32	41	4	18	7	.038	65	11, 22, 34
G07241	C-TQTS-3600	36	41	4	18	7	.038	65	11, 22, 34
Double lumen set									
G36224	C-DLCTS-2400	24	40	8	18	7	.038	80	12, 18, 26

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.





Thal-Quick Chest Tube Adapter

Used to provide separate sampling or an infusion port within Thal-Quick Chest Tubes.

Features and benefits

- The adapter is compatible with 8-36 Fr Thal-Quick Chest Tubes.

Order Number	Reference Part Number
G07467	C-TQF-28.0

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.



Fuhrman Pleural/Pneumopericardial Drainage Set

Used for evacuation of air from the pericardial sac or to drain air or fluid from the pleural space.





Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The catheter is available in a variety of sizes, ranging from 5-12 Fr, to cater to different patient anatomies.
- Sideports are positioned within the catheter pigtail and are designed to improve the drainage rate.
- Centimetre markings on the pigtail catheter allow for accurate placement.
- The radiopaque catheter material enhances x-ray visualisation.
- Pigtail catheters have a comparable efficacy to chest tubes in patients with pneumothorax.¹
- The catheter is compatible with the Cook Chest Drain Valve, which is designed to permit continued re-expansion of the lung, preventing inadvertent backflow of air into the pleural space.

1.Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma Acute Care Surg.* 2011;71(5):1104-1107.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G55713	C-PPD-500-WCE-IMH	5	15	6	18	5	.028	40	4.5
G55714	C-PPD-600-WCE-IMH	6	15	6	18	5	.028	40	5.5
G55723	C-PPD-850-IMH	8.5	15	6	18	7	.028	40	8
G55715	C-PPD-850-WCE-IMH	8.5	15	6	18	5	.028	40	8
G55716	C-PPD-1020-WCE-IMH	10.2	15	6	18	5	.028	40	9.5
G55727	C-PPD-1200-CHILDRENS-082694-IMH	12	15	6	18	7	.028	40	11

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.





Wood Pleural/Pneumopericardial Drainage Set

Used for evacuation of air from the pericardial sac or drainage of air or fluid from the pleural space.







Features and benefits

- The catheter/needle assembly facilitates minimally invasive catheter introduction via the catheter-over-needle technique.
- Sideports are positioned within the catheter pigtail and are designed to improve the drainage rate.
- The radiopaque catheter material enhances x-ray visualisation.
- Pigtail catheters have a comparable efficacy to chest tubes in patients with pneumothorax.¹

1.Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma Acute Care Surg.* 2011;71(5):1104-1107.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm
G55718	C-PPDJ-500-WOOD-IMH	5	8	10	23	16
G55720	C-PPDJ-600-WOOD-IMH	6	8	10	22	16

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.



24 X 36 INCH DRAPE (NOT SHOWN)



Thoracentesis Set


Used to remove fluid from the pleural space.

Features and benefits


- Seldinger placement facilitates controlled, minimally invasive catheter introduction.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm
G03286	C-THS-850	8.5	22	4	18	10	.038	70



Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.



ONE-PART PERCUTANEOUS ENTRY NEEDLE WITH BASEPLATE



SAFE-T-J FIXED CORE WIRE GUIDE



TWO VINYL CONNECTING TUBES



Pericardiocentesis Set - Straight Catheter

Used to remove fluid from the pericardial sac.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The radiopaque catheter material enhances x-ray visualisation.
- The stainless steel hub of the entry needle set permits constant electrocardiographic monitoring via alligator clip cable.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm
G03282	C-PCS-850	8.5	22	4	18	15	.038	70


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




PERICARDIOCENTESIS CATHETER




SAFE-T-J FIXED CORE WIRE GUIDE




CARDIOCENTESIS NEEDLE



VINYL CONNECTING TUBE



VINYL CONNECTING TUBE WITH THREE-WAY STOPCOCK



ALLIGATOR CLIP CABLE



Pericardiocentesis Set - Pigtail Catheter





Used to remove fluid from the pericardial sac.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The radiopaque catheter material enhances x-ray visualisation.
- Sideports are positioned within the catheter pigtail and are designed to improve the drainage rate.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G43558	C-PCS-PHH	7.1	30	6	18	7	.038	60	6, 7

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.





Pericardiocentesis Set - Pigtail Catheter


Used to remove fluid from the pericardial sac.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The radiopaque catheter material enhances x-ray visualisation.
- Sideports are positioned within the catheter pigtail and are designed to improve the drainage rate.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr	Drainage Bag mL
G10643	C-PCS-700-TORONTO-042997	7	22	4	18	9	.038	70	8	1400


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




PERICARDIOCENTESIS CATHETER



SAFE-T-J FIXED CORE WIRE GUIDE




ONE-PART PERCUTANEOUS ENTRY NEEDLE



DILATOR



NEEDLE HOLDER CUP



DRAINAGE BAG

NOT PICTURED - CSR WRAP



Lock Pericardiocentesis Catheter Set - Pigtail Catheter


Used to remove fluid from the pericardial sac.

Features and benefits


- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The 40 cm long catheter is designed to cater to larger patient anatomies.
- Sideports are positioned on both sides of the catheter and are designed to improve the drainage rate.
- The radiopaque catheter material enhances x-ray visualisation.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G02433	C-PCS-830-LOCK	8.3	40	6	18	7, 12	.038	70	8


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.




LOCK PERICARDIOCENTESIS CATHETER



DILATOR



SAFE-T-J PTFE-COATED FIXED CORE WIRE GUIDE



TWO ONE-PART PERCUTANEOUS ENTRY NEEDLES



PLASTIC MULTIPURPOSE TUBING ADAPTER



Teitel Pericardiocentesis Set


Used to remove fluid from the pericardial sac.

Features and benefits


- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The stainless steel hub of the entry needle permits constant electrocardiographic monitoring via alligator clip cable.
- The radiopaque catheter material enhances x-ray visualisation.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G05251	C-PCS-500-TTL	5	30	6	19	6	.028	60	6

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.



TEITEL PERICARDIOCENTESIS CATHETER




ONE-PART PERCUTANEOUS ENTRY NEEDLE




SAFE-T-J PTFE-COATED FIXED CORE WIRE GUIDE



DILATOR



TWO VINYL CONNECTING TUBES



ALLIGATOR CLIP CABLE

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



CC, URO-A4