

Melker Cuffed Emergency Cricothyrotomy Catheter Set Seldinger

Intended to establish emergency airway access when endotracheal intubation cannot be performed.



Features and benefits


- Airway access is achieved by using the Seldinger technique via the cricothyroid membrane.
- The airway catheter is cuffed to protect and control the airway after it is placed.

Components

- Melker Cuffed Emergency Cricothyrotomy Catheter
- Introducer needle
- PTFE catheter introducer needle
- 6 mL syringe
- #15 mini scalpel
- Amplatz Extra-Stiff Support Wire Guide
- Curved dilator
- Tracheostomy tape
- Needle holder cup

Order Number	Reference Part Number	Catheter Length cm	Catheter ID mm	Catheter OD mm	Wire Diameter inch	Wire Length cm	Needle gage	Needle Length cm	Quantity
G29432	C-TCCSB-500	9	5	7.2	.038	40	18	7	1


Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.



Melker Emergency Cricothyrotomy Catheter Set

Seldinger

Intended to establish emergency airway access when endotracheal intubation cannot be performed.



Features and benefits

- Airway access is achieved by using the Seldinger technique via the cricothyroid membrane.
- The uncuffed airway catheter comes in three different sizes to address a variety of patient anatomies.

Components

- Melker Emergency Cricothyrotomy Catheter
- Introducer needle
- PTFE catheter introducer needle
- 6 mL syringe
- #15 mini scalpel
- Amplatz Extra-Stiff Support Wire Guide
- Curved dilator
- Tracheostomy tape

Order Number	Reference Part Number	Catheter Length cm	Catheter ID mm	Catheter OD mm	Wire Diameter inch	Wire Length cm	Needle gage	Needle Length cm	Quantity
G06245	C-TCCS-350	3.8	3.5	5.1	.038	30	18	4	1
G06246	C-TCCS-400	4.2	4	5.9	.038	30	18	4	1
G04652	C-TCCS-600	7.5	6	8.8	.038	40	18	7	1

With tracheostomy tube plug

G07546	C-TCCS-400-TT	7.5	4	5.9	.038	40	18	7	1
--------	---------------	-----	---	-----	------	----	----	---	---

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.