

**REDPINE**

Endo Leads Better Life

# Single-Use Video Flexible Ureterorenoscope

RP-U-C0305

**7.5 Fr**



Excellent Irrigation

High Performance with Bendable Sheath



# Single-Use Video Flexible Ureterorenoscope



## Superior Optics Design

Unique transparent distal CMOS and LED design brings higher illumination



## Excellent Irrigation

Slim shaft with 3.6Fr working channel





## High performance

In the bendable sheath, bending section still performs well



## Smooth Maneuverability

Lightweight handle ergonomically designed makes better maneuverability



Model	RP-U-C0305
<b>Optical System</b>	
Field of view	120°
Depth of field	5mm-50mm
Direction of view	0°
<b>Insertion Portion</b>	
Outer diameter	Φ2.5mm(7.5Fr)
Inner diameter	3.6Fr
Working length	670mm (26.4in)
Bending angle	275°up, 275°down
<b>Other Specifications</b>	
Bending angle control	Angle lock

## Endoscopic Video Image Processor



RP-IPD-V800



### Full-HD Display

10.1-inch touchscreen provides synchronized display connections at ease



### Patient File Management System

Video and image during operation can be effortlessly collected and managed for data analysis



### Excellent Compatibility

Compatible with multiple REDPINE single-use endoscopes



### Multi-Functional Image Processor

Can switch between multiple endoscopic modes, personalize different visual frames and adjust image color and direction

**REDPINE**

Endo Leads Better Life

Guangzhou Red Pine Medical instrument Co., Ltd.

Tel: +86 20 80927258

Email: info@gzredpine.com

http://en.redpinemed.com

REDPINE\_Medical



Copyright © 2015 REDPINE

RP-U-C0305-202408-V.1