

# ENT patient chairs & working stools

medseat type





Since 1983, **Heinemann Medizintechnik** is one of the leading equipment providers in **ENT medical technology**.

The company is not only a **manufacturer** of high-quality ENT treatment units, but also a **competent** partner for practice planning as well as your sales **partner** for complete ENT equipment in **hospitals** and private practices.

# Patient chairs

Medseat, with its numerous variants, is the most popular treatment chair in ENT practices worldwide. It has won the trust of many practitioners through its ability to provide high flexibility and variability in a robust form that is easy to use. In addition, we have redesigned it over and over again throughout its product life. Today, in view of the growing scope of treatment in practices, the need for ergonomics for doctor and patient is also growing. It is about much more than just sitting upright, lying comfortably or maintaining a healthy posture for medical diagnosis and treatment. The focus is much more on the features of the treatment chair with an ergonomically movable fold-out step and flexible headrests. medseat is a compact all-rounder. It takes up little space and is suitable for examinations, advice and treatment, including surgical procedures.

In ENT practices, it is an integral part of many practice solutions. But it is also used in dentist, orthodontic or oral surgery practices as well as in plastic surgery.

More than ever, the new generation of the medseat accommodates your individuality in terms of treatment: Because now you can configure your treatment chair yourself combined with the head and footrest of your choice. For example, you can decide, tailor-made for your needs, whether you want a synchronized leg section with a guided Ergo fold-out step and the Ergo headrest, whether you need the Vario double joint with a secure hold, or whether you prefer to work with a foot switch or hand switch in the backrest.

The optional Quick Assist quick coupling for the head restraint, half or full spherical caps is a very ergonomic and efficient solution for frequent changes in reclining and positioning the patient. These variants only outline a selection of your options, in addition to the choice of colours and covers and a wide range of accessories. medseat is a classic that always stays current with new ideas.

















<b>Medseat Models</b>	<b>REF</b>
Medseat 1.1 ENT patient chair, hydraulic (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4612050
Medseat 1.2 ENT patient chair, hydraulic with synchronized legrest and Ergo headrest	553-4632000
Medseat 2.1 ENT patient chair, electromechanical (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4616150
Medseat 2.2 ENT patient chair, electromechanical with synchronized legrest and Ergo headrest	553-4636100
Medseat 2.3 ENT patient chair, electromechanical with hand control in the backrest and Ergo headrest	553-4616108
Medseat 2.4 ENT patient chair, electromechanical with hand control in the backrest, synchronized legrest and Ergo headrest	553-4636108
Medseat 3.1 ENT patient chair, electromechanical (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4716150
Medseat 3.2 ENT patient chair with, electromechanical synchronized legrest and Ergo headrest	553-4736100
Medseat 3.3 ENT patient chair, electromechanical with hand control in the backrest (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4716158
Medseat 3.4 ENT patient chair, electromechanical with hand control in the backrest, synchronized legrest and Ergo headrest	553-4736108
Medseat 3.5 ENT patient chair, electromechanical with hand control in the backrest, legrest and Ergo headrest	553-4736118
Medseat 3.6 ENT patient chair, electromechanical with hand control in the backrest, synchronized legrest with integrated Ergo fold-out step (headrest not included, configurable at buyer's option - see accessories)	553-4736158
Medseat 2.5 ENT patient chair, electromechanical (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4616250
Medseat 2.6 ENT patient chair, electromechanical with synchronized legrest and Ergo headrest	553-4636200
Medseat 2.7 ENT patient chair, electromechanical with hand control in the backrest and Ergo headrest	553-4616208
Medseat 2.8 ENT patient chair, electromechanical with hand control in the backrest, synchronized legrest and Ergo headrest	553-4636208
Medseat 3.7 ENT patient chair, electromechanical (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4716250
Medseat 3.8 ENT patient chair, electromechanical with synchronized legrest and Ergo headrest	553-4736200
Medseat 3.9 ENT patient chair, electromechanical with hand control in the backrest (headrest and legrest not included, configurable at buyer's option - see accessories)	553-4716258
Medseat 3.10 ENT patient chair, electromechanical with hand control in the backrest, synchronized legrest and Ergo headrest	553-4736208
Medseat 3.11 ENT patient chair, electromechanical with hand control in the backrest, legrest and Ergo headrest	553-4736218
Medseat 4.1 ENT patient chair, electromechanical with hand control in the backrest, synchronized legrest with integrated Ergo fold-out step (headrest not included, configurable at buyer's option - see accessories)	553-4736258





Model	Height adjustment		Upper part		Footrest		Back rest	Base
	hydraulic	electrical	manual	electrical	synchronous	electrical (separately)	hand control	mobile
Medseat 1.1	X		X					
Medseat 1.2	X		X		X			
Medseat 2.1		X	X					
Medseat 2.2		X	X		X			
Medseat 2.3		X	X				X	
Medseat 2.4		X	X		X		X	
Medseat 3.1		X		X				
Medseat 3.2		X		X	X			
Medseat 3.3		X		X			X	
Medseat 3.4		X		X	X		X	
Medseat 3.5		X		X		X	X	
Medseat 3.6		X		X	X		X	
Medseat 2.5		X	X					X
Medseat 2.6		X	X		X			X
Medseat 2.7		X	X				X	X
Medseat 2.8		X	X		X		X	X
Medseat 3.7		X		X				X
Medseat 3.8		X		X	X			X
Medseat 3.9		X		X			X	X
Medseat 3.10		X		X	X		X	X
Medseat 3.11		X		X		X	X	X
Medseat 4.1		X		X	X		X	X

Options for Medseat Models	REF	
Additional charge Medseat in premium fabric (s. page 10)	553-0000461	
Additional charge Medseat in special colours	553-0000008	
Electrically conductive cover (black)	553-000EL11	
Ergo headrest	553-1461100	
Ergo headrest with Vario double joint	553-1465000	
Quick assist, quick coupling for various head bolsters	553-1465005	
Replaceable upholstery for Quick assist	553-1465011	
Full spherical cap for Quick assist	553-1465012	
Half spherical cap for Quick assist	553-1465013	
Additional charge swivable seat cushion	553-DSHNO	
Child's seat inlay	553-9120040	
Transparent cover leg part (retrofit not possible)	553-1414630	
Synchronized legrest with step	553-1335000	
synchronized legrest with integrated Ergo fold-out step	553-1335500	
Legrest fixed with step	553-1330904	
Legrest fixed with fold-out step	553-1337100	
Footrest bracket	553-1330850	

# Working chairs

---

The variety of the GREINER work stools programme is as versatile as the working life in the medical field. Well known for their „Made in Germany“ experience, the stools have been supporting the work of medical professionals 24/7 around the world for years. The range of applications extends from examination and treatment up to the operating room. Both because of their strengths in anatomy, ergonomics, and functionality and because of accommodating GREINER medical chairs for specialist applications. In other words: day after day, GREINER helps to make the everyday work of healthcare better, easier, and healthier. Good seating in the best quality which facilitates processes while strengthening the back.



With innovative features, GREINER has always managed to show that there is a difference in work stools. This includes the fulfilment of hygienic and technical requirements for treatment as well as the varied selection of heights for sitting and upholstery, height adjustment with a 360° manual ring or foot release of the gas spring as well as the choice of base part and castors. Work stools and chairs of high-class quality that retain a long working life are produced at the in-house facilities including a foaming and sewing department. The cover in premium material is extremely comfortable, which noticeably improves the sitting comfort. Of course, all work stools meet the clinical requirements for stability, safety and hygiene, while sustainability in production and permanent spare parts supply are guaranteed.






<b>Work stool models</b>	<b>REF</b>
Work stool 1.1, 40-52 cm, ring release, seat ø 36 cm	553-3008105
Work stool 1.1, 46-58 cm, ring release, seat ø 36 cm	553-3008205
Work stool 1.1, 54-73 cm, ring release, seat ø 36 cm	553-3008305
Work stool 1.2, 46-58 cm, ring release, seat ø 36 cm, with foot ring	553-3009005
Work stool 1.2, 54-73 cm, ring release, seat ø 36 cm, with foot ring	553-3009105
Work stool 2.1, 40-52 cm, ring release, seat ø 40 cm	553-4108105
Work stool 2.1, 46-58 cm, ring release, seat ø 40 cm	553-4108205
Work stool 2.1, 54-73 cm, ring release, seat ø 40 cm	553-4108305
Work stool 2.2, 46-58 cm, ring release, seat ø 40 cm, with foot ring	553-4109005
Work stool 2.2, 54-73 cm, ring release, seat ø 40 cm, with foot ring	553-4109105
Work stool 3.1, 40-52 cm, ring release, anatomically shaped seat ø 40 cm	553-4128105
Work stool 3.1, 46-58 cm, ring release, anatomically shaped seat ø 40 cm	553-4128205
Work stool 3.1, 54-73 cm, ring release, anatomically shaped seat ø 40 cm	553-4128305
Work stool 3.2, 46-58 cm, ring release, anatomically shaped seat ø 40 cm, with foot ring	553-4129005
Work stool 3.2, 54-73 cm, ring release, anatomically shaped seat ø 40 cm, with foot ring	553-4129105
Work stool 4.1, 43-56 cm, foot release, seat ø 40 cm, with foot ring	553-4105205
Work stool 4.1, 51-71 cm, foot release, seat ø 40 cm, with foot ring	553-4105305
Work stool 4.2, 43-56 cm, foot release, anatomically shaped seat ø 40 cm, with foot ring	553-4125205
Work stool 4.2, 51-71 cm, foot release, anatomically shaped seat ø 40 cm, with foot ring	553-4125305



<b>Work chair models</b>	<b>REF</b>
Work chair 1.1, 40-52 cm, small backrest, ring release, seat ø 36 cm	553-3018105
Work chair 1.1, 46-58 cm, small backrest, ring release, seat ø 36 cm	553-3018205
Work chair 1.1, 54-73 cm, small backrest, ring release, seat ø 36 cm	553-3018305
Work chair 1.2, 46-58 cm, small backrest, ring release, seat ø 36 cm, with foot ring	553-3019005
Work chair 1.2, 54-73 cm, small backrest, ring release, seat ø 36 cm, with foot ring	553-3019105
Work chair 2.1, 40-52 cm, large backrest, ring release, seat ø 40 cm	553-3028105
Work chair 2.1, 46-58 cm, large backrest, ring release, seat ø 40 cm	553-3028205
Work chair 2.1, 54-73 cm, large backrest, ring release, seat ø 40 cm	553-3028305
Work chair 2.2, 46-58 cm, large backrest, ring release, seat ø 40 cm, with foot ring	553-3029005
Work chair 2.2, 54-73 cm, large backrest, ring release, seat ø 40 cm, with foot ring	553-3029105
Work chair 3.1, 43-53 cm, crescent-shaped backrest, manual ring release, seat ø 40 cm	553-4154105
Work chair 3.1, 47-59 cm, crescent-shaped backrest, manual ring release, seat ø 40 cm	553-4154205
Work chair 3.1, 54-73 cm, crescent-shaped backrest, manual ring release, seat ø 40 cm	553-4154305
Work chair 3.2, 47-59 cm, crescent-shaped backrest, manual ring release, seat ø 40 cm, with foot ring	553-4155005
Work chair 3.2, 54-73 cm, crescent-shaped backrest, manual ring release, seat ø 40 cm, with foot ring	553-4155105
Work chair 4.1, 40-52 cm, crescent-shaped backrest, manual ring release, anatomically shaped seat ø 40 cm	553-4184105
Work chair 4.1, 46-58 cm, crescent-shaped backrest, manual ring release, anatomically shaped seat ø 40 cm	553-4184205
Work chair 4.1, 54-73 cm, crescent-shaped backrest, manual ring release, anatomically shaped seat ø 40 cm	553-4184305
Work chair 4.2, 46-58 cm, crescent-shaped backrest, manual ring release, anatomically shaped seat ø 40 cm, with foot ring	553-4185005
Work chair 4.2, 54-73 cm, crescent-shaped backrest, manual ring release, anatomically shaped seat ø 40 cm, with foot ring	553-4185105
Work chair 4.3, 43-56 cm, crescent-shaped backrest, foot release, seat ø 40 cm, with foot ring	553-4155205
Work chair 4.3, 51-71 cm, crescent-shaped backrest, foot release, seat ø 40 cm, with foot ring	553-4155305
Work chair 4.4, 43-56 cm, crescent-shaped backrest, foot release, anatomically shaped seat ø 40 cm, with foot ring	553-4185205
Work chair 4.4, 51-71 cm, crescent-shaped backrest, foot release, anatomically shaped seat ø 40 cm, with foot ring	553-4185305
Work chair 5.1, 43-56 cm, blarge backrest, foot release, seat ø 40 cm, with foot ring	553-3045205
Work chair 5.1, 51-71 cm, large backrest, foot release, seat ø 40 cm, with foot ring	553-3045305



Options for work stool/ -chair models	REF	
additional charge work stool in premium fabric (see below)	553-0000462	
additional charge work chair in premium fabric (see below)	553-0000463	
electrically conductive cover (black) for work stools	553-001EL11	
electrically conductive cover (black) for work chairs	553-002EL11	
round plate (cover in stainless steel)	553-0000464	
Glider instead of castors	553-000RG01	
central arrest	553-0000465	

**PREMIUM Upholstery colours** (optional against surcharge):



P800 (white)



P200 (lightgrey)




P210 (anthracite)



P400 (grass)



P300 (mocca)



P310 (chocolate)



P500 (citron)



P600 (carmin)



P610 (fuchsia)



P700 (aqua)



P710 (navy)



P100 (black)


## STANDARD Upholstery colours:



C811 (white)



C450 (mint)



C422 (cyan)



C435 (petrol)



C213 (lightgrey)



C214 (stone)



C220 (anthracite)



C111 (black)



C830 (beige)



C350 (brown)



C330 (mocca)



C360 (chocolate)



C700 (sky)



C706 (atoll)



C711 (navy)



C715 (purple)



C500 (citron)



C530 (orange)



C637 (rosso)



C610 (fuchsia)

# Contact

G. Heinemann Medizintechnik GmbH

Leibnizstraße 13-15

24568 Kaltenkirchen

Germany

Phone: +49 (0)4191 - 95379-0

Fax: +49 (0)4191 - 9537937

Email: [info@heinemann-ent.de](mailto:info@heinemann-ent.de)

Web: [www.heinemann-ent.de](http://www.heinemann-ent.de)